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How to Use this Product

Guidance on use in the KAPA HyperCap Workflow can be found in the KAPA HyperChoice, KAPA HyperExplore, and KAPA HyperExome Instructions for Use.

	1	2	3	4	5	6	7	8	9	10	11	12
A	UDI-P 01	UDI-P 09	UDI-P 17	UDI-P 25	UDI-P 33	UDI-P 41	UDI-P 49	UDI-P 57	UDI-P 65	UDI-P 73	UDI-P 81	UDI-P 89
B	UDI-P 02	UDI-P 10	UDI-P 18	UDI-P 26	UDI-P 34	UDI-P 42	UDI-P 50	UDI-P 58	UDI-P 66	UDI-P 74	UDI-P 82	UDI-P 90
C	UDI-P 03	UDI-P 11	UDI-P 19	UDI-P 27	UDI-P 35	UDI-P 43	UDI-P 51	UDI-P 59	UDI-P 67	UDI-P 75	UDI-P 83	UDI-P 91
D	UDI-P 04	UDI-P 12	UDI-P 20	UDI-P 28	UDI-P 36	UDI-P 44	UDI-P 52	UDI-P 60	UDI-P 68	UDI-P 76	UDI-P 84	UDI-P 92
E	UDI-P 05	UDI-P 13	UDI-P 21	UDI-P 29	UDI-P 37	UDI-P 45	UDI-P 53	UDI-P 61	UDI-P 69	UDI-P 77	UDI-P 85	UDI-P 93
F	UDI-P 06	UDI-P 14	UDI-P 22	UDI-P 30	UDI-P 38	UDI-P 46	UDI-P 54	UDI-P 62	UDI-P 70	UDI-P 78	UDI-P 86	UDI-P 94
G	UDI-P 07	UDI-P 15	UDI-P 23	UDI-P 31	UDI-P 39	UDI-P 47	UDI-P 55	UDI-P 63	UDI-P 71	UDI-P 79	UDI-P 87	UDI-P 95
H	UDI-P 08	UDI-P 16	UDI-P 24	UDI-P 32	UDI-P 40	UDI-P 48	UDI-P 56	UDI-P 64	UDI-P 72	UDI-P 80	UDI-P 88	UDI-P 96

Fig. 1: Layout of the KAPA UDI Primer Mixes Plate.

Pooling/Multiplexing Guidelines

As a rule, choose primer mixes as highlighted by the boxes in Figure 2 to take advantage of the color-balanced indexes. This will prevent registration failure and laser color complexity issues during sequencing and de-multiplexing. Figure 2 details the recommended two-plexing combinations that are fully color-balanced.

See the detailed suggestions below:

- Pooling two samples (two-plex):
 - Figure 2 demonstrates the recommended two-plex combinations that are fully color-balanced based on the plate layout in Figure 1 (all combinations indicated by the gray boxes, e.g. A1 + B1, C1 + D1, etc.).
- Pooling three to eight samples (three to eight-plex):
 - Use any of the recommended two-plex combinations with any other index in the column (all combinations indicated by colored box e.g. three-plex: A1 + B1 + C1, four-plex: A1 + B1 + C1 + D1, etc.).
- Pooling nine or more samples:
 - Any number of additional libraries may be multiplexed with any of the color-balanced combinations listed below to obtain pools of any plexity.
 - It is recommended to use column groups of indexes (e.g. colored box 1, 2 etc.).


Fig. 2: Recommendations for two-plex (by pairs in grayscale boxes) and eight-plex (by columns in blue and green boxes) pooling. Note that the notched corner is on the bottom left of the plate. This directs the correct orientation of the plate with A1 positioned in the top left of the plate.

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Chapter 4. Amplify the Sample Library Using the KAPA UDI Primer Mixes

This chapter describes how to amplify the sample library using the KAPA UDI Primer Mixes in preparation for hybridization to the KAPA Target Enrichment Probes. This chapter requires the use of the components from the following kits:

Step 1. Prepare the Pre-Capture PCR Reaction

-  We recommend the inclusion of negative (water) and positive (previously amplified library) controls in the Pre-Capture PCR step.
-  For guidance on pre-capture and post-capture sample multiplexing, please refer to the KAPA UDI Primer Mixes Instructions for Use.
-  Make sure the KAPA HyperPure Beads are removed from storage to allow time for proper equilibration to room temperature. For best performance, store the beads protected from light when not in use.
-  Ensure to record the well position of the KAPA UDI Primer Mixes used for each sample.
 1. Retrieve and thaw the KAPA UDI Primer Mixes plate prepared in [Chapter 2 Step 3](#).
 2. Spin the plate at 280 x g for 30 seconds to collect the contents to the bottom of the wells.
 3. Peel off or pierce the foil seal for the appropriate number of wells needed.
 -  If piercing the foil seal, avoid cross contamination by using a new pipette tip for every well.
 4. Add 5 µL of a KAPA UDI Primer Mixes to each individual sample library.
 5. Add 25 µL of KAPA HiFi Hotstart ReadyMix to each combined sample library and KAPA UDI Primer Mixes.
 6. Mix thoroughly and perform a quick spin. Immediately proceed to amplification.

Step 2. Perform the Pre-Capture PCR Amplification

7. Place the sample in the thermocycler and amplify the sample library using the following Pre-Capture PCR program with the lid temperature set to +105°C:
 - Step 1: 45 seconds at +98°C
 - Step 2: 15 seconds at +98°C
 - Step 3: 30 seconds at +60°C
 - Step 4: 30 seconds at +72°C
 - Step 5: Go to Step 2, Variable (see table below for recommendation)
 - Step 6: 1 minute at +72°C
 - Step 7: Hold at +4°C

Library Preparation Kit	Go to Step 2:
KAPA HyperPrep Kit	7 times (8 total cycles)
KAPA HyperPlus Kit	5 times (6 total cycles)

8. Proceed immediately to the next step.

Step 3. Purify the Amplified Sample Library using KAPA HyperPure Beads

1. Add 70 µL of room temperature, thoroughly resuspended, KAPA HyperPure Beads to each amplified sample library.
2. Mix the amplified sample library and KAPA HyperPure Beads thoroughly and perform a quick spin.

It is important at this step to ensure that the solution is thoroughly mixed and appears homogenous. Insufficient mixing may compromise recovery.

3. Incubate the sample at room temperature for 5 minutes to allow the DNA to bind to the beads.
4. Place the sample on a magnet to collect the beads. Incubate until the liquid is clear.
5. Carefully remove and discard the supernatant.
6. Keeping the sample on the magnet, add 200 µL of freshly-prepared 80% ethanol.
7. Incubate the sample at room temperature for ≥30 seconds.
8. Carefully remove and discard the ethanol.
9. Keeping the sample on the magnet, add 200 µL of freshly-prepared 80% ethanol.
10. Incubate the sample at room temperature for ≥30 seconds.
11. Carefully remove and discard the ethanol. Remove residual ethanol without disturbing the beads.
12. Allow the beads to dry at room temperature, sufficiently for all of the ethanol to evaporate.

Over-drying the beads may result in dramatic yield loss.

13. Remove the sample from the magnet.
14. Thoroughly resuspend the beads in 32 µL of 10 mM Tris-HCl, pH 8.0 or PCR Grade water.
15. Incubate the sample at room temperature for 2 minutes to allow the DNA to elute off the beads.
16. Place the sample on a magnet to collect the beads. Incubate until the liquid is clear.
17. Transfer 30 µL of the eluate to a new tube/well.
18. Purified, amplified libraries can be stored at +2 to +8°C for 1-2 weeks or at -15 to -25°C for long term storage.

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KAPA RNA HyperPrep Kit with RiboErase (HMR)
Illumina® Platforms

Product Applications

The KAPA RNA HyperPrep Kit with RiboErase (HMR) is designed for both manual and automated NGS library construction from 25 ng – 1 µg of total RNA. The kit depletes both cytoplasmic (5S, 5.8S, 18S, and 28S), and mitochondrial (12S and 16S) rRNA species. The protocol is applicable to a wide range of RNA-Seq applications, including:

- gene expression analysis of high- and low-quality RNA samples (e.g., extracted from FFPE tissue);
- single nucleotide variation (SNV) discovery;
- splice junction and gene fusion identification; and
- characterization of both polyadenylated and non-polyadenylated RNAs, including noncoding and immature RNAs.

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