

Plec de prescripcions tècniques per a la contractació del Subministrament de reactius per als tests de detecció genètica per l'Institut Català d'Oncologia

1. Objecte de la contractació

L'objecte d'aquesta licitació és la contractació del subministrament de tests de detecció genètica per l'Institut Català d'Oncologia, en base a les quantitats que consten en el Informe de necessitats. Aquests tests són reactius i sondes de MLPA®, utilitzats per analitzar la possible variació en el nombre de còpies de nombrosos gens. Aquesta tècnica està considerada avui dia com el "gold standard" per a aquest tipus de determinacions genètiques. El MLPA és una tècnica que permet en una mateixa reacció múltiplex analitzar el nombre de còpies de diverses regions d'una manera semiquantitatiu, es basa en una primera reacció d'unió-ligació de sondes amb la zona homòloga d'interès i la posterior amplificació de les sondes lligades. Mitjançant un anàlisi de fragments posterior en un seqüenciador capil·lar i aprofitant la diferència de mida de cadascuna de les sondes, es poden identificar aberracions en el nombre de còpies. Per realitzar aquesta tècnica es necessiten uns reactius comuns (SALSA MLPA EK5 reagent kit) i altres (les sondes) que són específics per a cada gen d'interès i que es poden adquirir en formats de 25, 50 i 100 reaccions. S'adjunta llistat de reactius de MLPA com annex núm. 1 al present plec de prescripcions tècniques.

Les quantitats de consum indicades a l'informe de necessitats es refereix a un any previstes per a cada article i que han estat extretes a partir dels consums històrics en el nostre sistema.

Aquestes quantitats es podran modificar a l'alça o a la baixa en funció de les necessitats de cada centre, depenent de la seva activitat assistencial, i d'acord amb allò que preveu la legislació vigent, sense superar el volum màxim previst. Els contractes que derivin de l'expedient de licitació són de subministraments de tracte successiu, sense que hi hagi un compromís de compra mínima. En cas de variació de quantitats, les empreses adjudicatàries no tindran dret a percebre cap indemnització per la variació.

2. Òrgan de contractació i centres destinataris del subministrament

L'òrgan de contractació és el Consell d'Administració de l'Institut Català d'Oncologia.

3. Estructura del procediment

Els materials a licitar s'estructuren en codis individuals. Els licitadors només podran presentar una oferta per a cadascun dels articles al qual concorrin. No s'admeten variants.

Així doncs, una única referència podrà ser només presentada contra un article. No és obligatori presentar oferta a tots i cadascun dels lots, en cas d'existir divisió en lots a l'informe a l'informe de necessitats i memòria justificativa. El licitador presentarà oferta per a cada format de 25, 50, 100 reaccions i SALSA MLPA EK5 reagent kit (540001484).

4. Condicions generals del subministrament

1. Els lliuraments sempre seran a ports pagats per l'adjudicatari.
2. Els productes objecte d'aquest plec seran conservats en gel i neu carbònica.
3. Es subministrarà el material sota comanda de l'Institut Català d'Oncologia.
4. El termini de lliurament serà, com a màxim, de 10 dies a comptats des de la data de recepció del fax o e-mail de les comandes corresponents. L'incompliment d'aquests terminis comportarà les sancions que s'especifiquen al plec de clàusules administratives particulars, podent fins i tot, comportar la resolució anticipada del contracte a instància de l'òrgan de contractació. Excepcionalment, la Unitat de Compres podrà valorar i autoritzar la modificació del termini de lliurament fixat per als productes objecte d'aquest procediment.
5. La data de caducitat de tots els productes referenciats en aquesta licitació no podrà ser inferior a 12 mesos.

Els productes oferts hauran de complir la normativa espanyola i comunitària vigent en matèria de qualitat, etiquetatge i envasat. En el cas de tractar-se de productes sanitaris, l'empresa licitadora haurà de presentar de forma imprescindible, la documentació acreditativa de compliment de la legislació vigent i de que els productes disposen del marcat de conformitat CE. Així mateix, pel que fa als productes sanitaris i les dades proporcionades pel fabricant, l'empresa haurà de presentar una declaració conforme tots els productes van etiquetats i aporten la documentació acomplint el contingut previst a l'apartat 13 de l'Annex I del RD 1591/2009, de 16 d'octubre, pel qual es regulen els productes sanitaris. Així mateix, els proveïdors hauran de proporcionar les acreditacions de complir amb les següents certificacions UNE-EN ISO 9001.

D'altra banda, les empreses adjudicatàries dels materials detallats en els dos paràgrafs anteriors seran responsables del manteniment de les caducitats dels mateixos, i de la reposició dels materials caducats, sense cap cost econòmic per l'ICO.

6. El lloc de lliurament serà els magatzems logístics dels que disposa l'Institut Català d'Oncologia, dins de l'horari del mateix, de 08:30 a 14:00h. de dilluns a divendres. El lliurament es realitzarà, si s'escau, en Palé europeu retractilat, amb un màxim d'1.40 cm d'altura.

7. Els productes de Material Sanitari que es detallen en aquest plec de prescripcions tècniques hauran de presentar-se en els envasos habituals de consum, per a cada tipus de producte. En cap cas les unitats d'embalatge hauran de superar el pes de 25 kg.

L'empresa adjudicatària encarregada de realitzar el subministrament, formalitzarà un albarà de lliurament del subministrament en el qual especificarà:

- Número de Comanda
- Data de lliurament
- Quantitat, tipus de producte i referència del producte lliurat.

Valoració econòmica del subministrament amb partida d'IVA inclosa.

5. Requeriments addicionals obligatoris

5.1 Formació, si el centre destinatari ho requereix, les empreses adjudicatàries hauran de formar al personal per al correcte ús dels productes subministrats.

5.2 Manual d'usuari: el proveïdor haurà de lliurar un manual d'usuari, en aquells casos que sigui necessari, document d'instruccions o protocol d'utilització del material, preferiblement en català, i si no es possible en castellà.

6. Oferta tècnica

Les ofertes tècniques consistiran en la Fitxa tècnica del producte.

És obligatori presentar les fitxes de dades de seguretat dels productes quan procedeixi. En el supòsit que les empreses licitadores no presentin la documentació tècnica del producte, amb la suficient informació i la forma de presentació no sigui la sol·licitada, l'oferta tècnica no serà avaluada i per tant es rebutjarà.

No és necessari incloure documentació bibliogràfica dels productes oferts. En cas que l'òrgan de contractació ho consideri oportú es requerirà a les empreses.

6.2 Presentació de mostres

És necessari aportar 3 mostres del material ofert, si així s'indica al quadre de característiques. L'òrgan de contractació es reserva el dret de sol·licitar-les durant el període d'avaluació de les ofertes si ho considera oportú.

En cas que es requereixin, el licitador haurà de lliurar les mostres sol·licitades al centre que es determini en un termini màxim de 24 hores.

En el supòsit que les empreses licitadores no presentin les mostres sol·licitades, l'oferta presentada no serà avaluada i, per tant, serà desestimada.

Amb les mostres presentades, l'Institut Català d'Oncologia es reserva el dret de fer-ne l'ús que estimi convenient.

7. Oferta econòmica

7.1 Preus de licitació

Els preus unitaris de licitació són els preus que figuren en l'informe de necessitats per a cada article.

L'oferta econòmica s'haurà de presentar, obligatòriament, en el model que s'adjunta com a Annex 1 al present plec de Clàusules Administratives Particulars.

8. Adjudicació

Amb caràcter general, l'adjudicació s'efectuarà per article, és a dir, de forma individual, recaient l'adjudicació de cadascun dels articles en un únic codi SAP i en un únic proveïdor.

9. Condicions logístiques i de subministrament

9.1 Condicions de subministrament

Les empreses adjudicatàries es comprometen a subministrar els materials adjudicats amb una caducitat mínima de 12 mesos. Les empreses adjudicatàries retiraran els materials defectuosos i els substituirà per altres de nous.

9.2 Garantia de subministrament

Les empreses adjudicatàries es comprometen a mantenir les existències que garanteixin una adequada continuïtat en l'abastiment a tots i cadascun dels centres destinataris. En el cas d'impossibilitat del subministrament:

- a) Hauran de comunicar-ho a la Unitat d'Aprovisionaments de l'Institut Català d'Oncologia.
- b) Facilitaran una alternativa pròpia amb l'antelació suficient, per tal de comprovar la seva idoneïtat i ser aprovada per part dels centres destinataris.
- c) En el supòsit que l'alternativa presentada no sigui equivalent al producte adjudicat inicialment, la Unitat d'Aprovisionaments prendrà les mesures oportunes per tal de garantir l'abastament als centres destinataris. Les empreses adjudicatàries es comprometen a garantir el subministrament durant la vigència del contracte d'acord amb les quantitats que s'estableixen en funció del consum dels centres. Les empreses adjudicatàries hauran de disposar d'un telèfon i/o correu electrònic específic de contacte per tal de poder tenir assegurat el subministrament o per qualsevol consulta que pogués sorgir. Així mateix, les empreses adjudicatàries es comprometen a assegurar els terminis de lliurament que s'indiquen en el present plec.

Jesús del Valle

Unitat de Diagnòstic Molecular

Institut Català d'Oncologia

L'Hospitalet de Llobregat, 20 de febrer de 2025

Llistat de reactius de MLPA (amb caràcter indicatiu):

Referència Descripció

ME001-025R	SALSA MLPA Probemix ME001 Tumour suppressor mix 1 - 25 reactions
ME001-050R	SALSA MLPA Probemix ME001 Tumour suppressor mix 1 - 50 reactions
ME001-100R	SALSA MLPA Probemix ME001 Tumour suppressor mix 1 - 100 reactions
ME002-025R	SALSA MLPA Probemix ME002 Tumour suppressor mix 2 - 25 reactions
ME002-050R	SALSA MLPA Probemix ME002 Tumour suppressor mix 2 - 50 reactions
ME002-100R	SALSA MLPA Probemix ME002 Tumour suppressor mix 2 - 100 reactions
ME011-025R	SALSA MLPA Probemix ME011 Mismatch Repair Genes - 25 reactions
ME011-050R	SALSA MLPA Probemix ME011 Mismatch Repair Genes - 50 reactions
ME011-100R	SALSA MLPA Probemix ME011 Mismatch Repair Genes - 100 reactions
ME012-025R	SALSA MLPA Probemix ME012 MGMT-IDH1-IDH2 - 25 reactions
ME012-050R	SALSA MLPA Probemix ME012 MGMT-IDH1-IDH2 - 50 reactions
ME012-100R	SALSA MLPA Probemix ME012 MGMT-IDH1-IDH2 - 100 reactions
ME024-025R	SALSA MLPA Probemix ME024 9p21 CDKN2A/2B region - 25 reactions
ME024-050R	SALSA MLPA Probemix ME024 9p21 CDKN2A/2B region - 50 reactions
ME024-100R	SALSA MLPA Probemix ME024 9p21 CDKN2A/2B region - 100 reactions
ME028-025R	SALSA MLPA Probemix ME028 PWS/AS - 25 reactions
ME028-050R	SALSA MLPA Probemix ME028 PWS/AS - 50 reactions
ME028-100R	SALSA MLPA Probemix ME028 PWS/AS - 100 reactions
ME029-025R	SALSA MLPA Probemix ME029 FMR1/AFF2 - 25 reactions
ME029-050R	SALSA MLPA Probemix ME029 FMR1/AFF2 - 50 reactions
ME029-100R	SALSA MLPA Probemix ME029 FMR1/AFF2 - 100 reactions
ME030-025R	SALSA MLPA Probemix ME030 BWS/RSS - 25 reactions
ME030-050R	SALSA MLPA Probemix ME030 BWS/RSS - 50 reactions
ME030-100R	SALSA MLPA Probemix ME030 BWS/RSS - 100 reactions
ME031-025R	SALSA MLPA Probemix ME031 GNAS - 25 reactions
ME031-050R	SALSA MLPA Probemix ME031 GNAS - 50 reactions
ME031-100R	SALSA MLPA Probemix ME031 GNAS - 100 reactions
ME032-025R	SALSA MLPA Probemix ME032 UPD7-UPD14 - 25 reactions
ME032-050R	SALSA MLPA Probemix ME032 UPD7-UPD14 - 50 reactions
ME032-100R	SALSA MLPA Probemix ME032 UPD7-UPD14 - 100 reactions
ME033-025R	SALSA MLPA Probemix ME033 TNDM - 25 reactions
ME033-050R	SALSA MLPA Probemix ME033 TNDM - 50 reactions
ME033-100R	SALSA MLPA Probemix ME033 TNDM - 100 reactions
ME034-025R	SALSA MLPA Probemix ME034 Multi-locus Imprinting - 25 reactions
ME034-050R	SALSA MLPA Probemix ME034 Multi-locus Imprinting - 50 reactions
ME034-100R	SALSA MLPA Probemix ME034 Multi-locus Imprinting - 100 reactions
ME042-025R	SALSA MLPA Probemix ME042 CIMP - 25 reactions
ME042-050R	SALSA MLPA Probemix ME042 CIMP - 50 reactions
ME042-100R	SALSA MLPA Probemix ME042 CIMP - 100 reactions
P002-025R	SALSA MLPA Probemix P002 BRCA1 - 25 reactions
P002-050R	SALSA MLPA Probemix P002 BRCA1 - 50 reactions
P002-100R	SALSA MLPA Probemix P002 BRCA1 - 100 reactions
P003-025R	SALSA MLPA Probemix P003 MLH1/MSH2 - 25 reactions
P003-050R	SALSA MLPA Probemix P003 MLH1/MSH2 - 50 reactions
P003-100R	SALSA MLPA Probemix P003 MLH1/MSH2 - 100 reactions
P008-025R	SALSA MLPA Probemix P008 PMS2 - 25 reactions
P008-050R	SALSA MLPA Probemix P008 PMS2 - 50 reactions
P008-100R	SALSA MLPA Probemix P008 PMS2 - 100 reactions
P010-025R	SALSA MLPA Probemix P010 POLG - 25 reactions
P010-050R	SALSA MLPA Probemix P010 POLG - 50 reactions

P010-100R	SALSA MLPA Probemix P010 POLG - 100 reactions
P011-025R	SALSA MLPA Probemix P011 VWF mix 1 - 25 reactions
P011-050R	SALSA MLPA Probemix P011 VWF mix 1 - 50 reactions
P011-100R	SALSA MLPA Probemix P011 VWF mix 1 - 100 reactions
P012-025R	SALSA MLPA Probemix P012 VWF mix 2 - 25 reactions
P012-050R	SALSA MLPA Probemix P012 VWF mix 2 - 50 reactions
P012-100R	SALSA MLPA Probemix P012 VWF mix 2 - 100 reactions
P013-025R	SALSA MLPA Probemix P013 ATRX - 25 reactions
P013-050R	SALSA MLPA Probemix P013 ATRX - 50 reactions
P013-100R	SALSA MLPA Probemix P013 ATRX - 100 reactions
P015-025R	SALSA MLPA Probemix P015 MECP2 - 25 reactions
P015-050R	SALSA MLPA Probemix P015 MECP2 - 50 reactions
P015-100R	SALSA MLPA Probemix P015 MECP2 - 100 reactions
P016-025R	SALSA MLPA Probemix P016 VHL - 25 reactions
P016-050R	SALSA MLPA Probemix P016 VHL - 50 reactions
P016-100R	SALSA MLPA Probemix P016 VHL - 100 reactions
P017-025R	SALSA MLPA Probemix P017 MEN1 - 25 reactions
P017-050R	SALSA MLPA Probemix P017 MEN1 - 50 reactions
P017-100R	SALSA MLPA Probemix P017 MEN1 - 100 reactions
P018-025R	SALSA MLPA Probemix P018 SHOX - 25 reactions
P018-050R	SALSA MLPA Probemix P018 SHOX - 50 reactions
P018-100R	SALSA MLPA Probemix P018 SHOX - 100 reactions
P021-025R	SALSA MLPA Probemix P021 SMA - 25 reactions
P021-050R	SALSA MLPA Probemix P021 SMA - 50 reactions
P021-100R	SALSA MLPA Probemix P021 SMA - 100 reactions
P022-025R	SALSA MLPA Probemix P022 PLP1 - 25 reactions
P022-050R	SALSA MLPA Probemix P022 PLP1 - 50 reactions
P022-100R	SALSA MLPA Probemix P022 PLP1 - 100 reactions
P025-025R	SALSA MLPA Probemix P025 Canavan - 25 reactions
P025-050R	SALSA MLPA Probemix P025 Canavan - 50 reactions
P025-100R	SALSA MLPA Probemix P025 Canavan - 100 reactions
P026-025R	SALSA MLPA Probemix P026 Sotos - 25 reactions
P026-050R	SALSA MLPA Probemix P026 Sotos - 50 reactions
P026-100R	SALSA MLPA Probemix P026 Sotos - 100 reactions
P027-025R	SALSA MLPA Probemix P027 Uveal melanoma - 25 reactions
P027-050R	SALSA MLPA Probemix P027 Uveal melanoma - 50 reactions
P027-100R	SALSA MLPA Probemix P027 Uveal melanoma - 100 reactions
P028-025R	SALSA MLPA Probemix P028 FHL - 25 reactions
P028-050R	SALSA MLPA Probemix P028 FHL - 50 reactions
P028-100R	SALSA MLPA Probemix P028 FHL - 100 reactions
P029-025R	SALSA MLPA Probemix P029 WBS - 25 reactions
P029-050R	SALSA MLPA Probemix P029 WBS - 50 reactions
P029-100R	SALSA MLPA Probemix P029 WBS - 100 reactions
P031-025R	SALSA MLPA Probemix P031 FANCA mix 1 - 25 reactions
P031-050R	SALSA MLPA Probemix P031 FANCA mix 1 - 50 reactions
P031-100R	SALSA MLPA Probemix P031 FANCA mix 1 - 100 reactions
P032-025R	SALSA MLPA Probemix P032 FANCA mix 2 - 25 reactions
P032-050R	SALSA MLPA Probemix P032 FANCA mix 2 - 50 reactions
P032-100R	SALSA MLPA Probemix P032 FANCA mix 2 - 100 reactions
P033-025R	SALSA MLPA Probemix P033 CMT1 - 25 reactions
P033-050R	SALSA MLPA Probemix P033 CMT1 - 50 reactions
P033-100R	SALSA MLPA Probemix P033 CMT1 - 100 reactions
P034-025R	SALSA MLPA Probemix P034 DMD-1 - 25 reactions
P034-050R	SALSA MLPA Probemix P034 DMD-1 - 50 reactions

P034-100R	SALSA MLPA Probemix P034 DMD-1 - 100 reactions
P035-025R	SALSA MLPA Probemix P035 DMD-2 - 25 reactions
P035-050R	SALSA MLPA Probemix P035 DMD-2 - 50 reactions
P035-100R	SALSA MLPA Probemix P035 DMD-2 - 100 reactions
P036-025R	SALSA MLPA Probemix P036 Subtelomeres Mix 1 - 25 reactions
P036-050R	SALSA MLPA Probemix P036 Subtelomeres Mix 1 - 50 reactions
P036-100R	SALSA MLPA Probemix P036 Subtelomeres Mix 1 - 100 reactions
P037-025R	SALSA MLPA Probemix P037 CLL-1 - 25 reactions
P037-050R	SALSA MLPA Probemix P037 CLL-1 - 50 reactions
P037-100R	SALSA MLPA Probemix P037 CLL-1 - 100 reactions
P038-025R	SALSA MLPA Probemix P038 CLL-2 - 25 reactions
P038-050R	SALSA MLPA Probemix P038 CLL-2 - 50 reactions
P038-100R	SALSA MLPA Probemix P038 CLL-2 - 100 reactions
P040-025R	SALSA MLPA Probemix P040 CLL - 25 reactions
P040-050R	SALSA MLPA Probemix P040 CLL - 50 reactions
P040-100R	SALSA MLPA Probemix P040 CLL - 100 reactions
P041-025R	SALSA MLPA Probemix P041 ATM-1 - 25 reactions
P041-050R	SALSA MLPA Probemix P041 ATM-1 - 50 reactions
P041-100R	SALSA MLPA Probemix P041 ATM-1 - 100 reactions
P042-025R	SALSA MLPA Probemix P042 ATM-2 - 25 reactions
P042-050R	SALSA MLPA Probemix P042 ATM-2 - 50 reactions
P042-100R	SALSA MLPA Probemix P042 ATM-2 - 100 reactions
P043-025R	SALSA MLPA Probemix P043 APC - 25 reactions
P043-050R	SALSA MLPA Probemix P043 APC - 50 reactions
P043-100R	SALSA MLPA Probemix P043 APC - 100 reactions
P044-025R	SALSA MLPA Probemix P044 NF2 - 25 reactions
P044-050R	SALSA MLPA Probemix P044 NF2 - 50 reactions
P044-100R	SALSA MLPA Probemix P044 NF2 - 100 reactions
P045-025R	SALSA MLPA Probemix P045 BRCA2/CHEK2 - 25 reactions
P045-050R	SALSA MLPA Probemix P045 BRCA2/CHEK2 - 50 reactions
P045-100R	SALSA MLPA Probemix P045 BRCA2/CHEK2 - 100 reactions
P046-025R	SALSA MLPA Probemix P046 TSC2 - 25 reactions
P046-050R	SALSA MLPA Probemix P046 TSC2 - 50 reactions
P046-100R	SALSA MLPA Probemix P046 TSC2 - 100 reactions
P047-025R	SALSA MLPA Probemix P047 RB1 - 25 reactions
P047-050R	SALSA MLPA Probemix P047 RB1 - 50 reactions
P047-100R	SALSA MLPA Probemix P047 RB1 - 100 reactions
P048-025R	SALSA MLPA Probemix P048 LMNA/MYOT/ZMPSTE24 - 25 reactions
P048-050R	SALSA MLPA Probemix P048 LMNA/MYOT/ZMPSTE24 - 50 reactions
P048-100R	SALSA MLPA Probemix P048 LMNA/MYOT/ZMPSTE24 - 100 reactions
P049-025R	SALSA MLPA Probemix P049 SLC6A8 - ABCD1 - 25 reactions
P049-050R	SALSA MLPA Probemix P049 SLC6A8 - ABCD1 - 50 reactions
P049-100R	SALSA MLPA Probemix P049 SLC6A8 - ABCD1 - 100 reactions
P050-025R	SALSA MLPA Probemix P050 CAH - 25 reactions
P050-050R	SALSA MLPA Probemix P050 CAH - 50 reactions
P050-100R	SALSA MLPA Probemix P050 CAH - 100 reactions
P051-025R	SALSA MLPA Probemix P051 Parkinson mix 1 - 25 reactions
P051-050R	SALSA MLPA Probemix P051 Parkinson mix 1 - 50 reactions
P051-100R	SALSA MLPA Probemix P051 Parkinson mix 1 - 100 reactions
P052-025R	SALSA MLPA Probemix P052 Parkinson mix 2 - 25 reactions
P052-050R	SALSA MLPA Probemix P052 Parkinson mix 2 - 50 reactions
P052-100R	SALSA MLPA Probemix P052 Parkinson mix 2 - 100 reactions
P054-025R	SALSA MLPA Probemix P054 FOXL2-TWIST1 - 25 reactions
P054-050R	SALSA MLPA Probemix P054 FOXL2-TWIST1 - 50 reactions

P054-100R	SALSA MLPA Probemix P054 FOXL2-TWIST1 - 100 reactions
P055-025R	SALSA MLPA Probemix P055 PAH - 25 reactions
P055-050R	SALSA MLPA Probemix P055 PAH - 50 reactions
P055-100R	SALSA MLPA Probemix P055 PAH - 100 reactions
P056-025R	SALSA MLPA Probemix P056 TP53 - 25 reactions
P056-050R	SALSA MLPA Probemix P056 TP53 - 50 reactions
P056-100R	SALSA MLPA Probemix P056 TP53 - 100 reactions
P057-025R	SALSA MLPA Probemix P057 FANCD2-PALB2 - 25 reactions
P057-050R	SALSA MLPA Probemix P057 FANCD2-PALB2 - 50 reactions
P057-100R	SALSA MLPA Probemix P057 FANCD2-PALB2 - 100 reactions
P058-025R	SALSA MLPA Probemix P058 IGHMBP2 - 25 reactions
P058-050R	SALSA MLPA Probemix P058 IGHMBP2 - 50 reactions
P058-100R	SALSA MLPA Probemix P058 IGHMBP2 - 100 reactions
P059-025R	SALSA MLPA Probemix P059 Dystonia - 25 reactions
P059-050R	SALSA MLPA Probemix P059 Dystonia - 50 reactions
P059-100R	SALSA MLPA Probemix P059 Dystonia - 100 reactions
P060-025R	SALSA MLPA Probemix P060 SMA Carrier - 25 reactions
P060-050R	SALSA MLPA Probemix P060 SMA Carrier - 50 reactions
P060-100R	SALSA MLPA Probemix P060 SMA Carrier - 100 reactions
P061-025R	SALSA MLPA Probemix P061 Lissencephaly - 25 reactions
P061-050R	SALSA MLPA Probemix P061 Lissencephaly - 50 reactions
P061-100R	SALSA MLPA Probemix P061 Lissencephaly - 100 reactions
P062-025R	SALSA MLPA Probemix P062 LDLR - 25 reactions
P062-050R	SALSA MLPA Probemix P062 LDLR - 50 reactions
P062-100R	SALSA MLPA Probemix P062 LDLR - 100 reactions
P064-025R	SALSA MLPA Probemix P064 Microdeletion Syndromes 1B - 25 reactions
P064-050R	SALSA MLPA Probemix P064 Microdeletion Syndromes 1B - 50 reactions
P064-100R	SALSA MLPA Probemix P064 Microdeletion Syndromes 1B - 100 reactions
P065-025R	SALSA MLPA Probemix P065 Marfan Syndrome-1 - 25 reactions
P065-050R	SALSA MLPA Probemix P065 Marfan Syndrome-1 - 50 reactions
P065-100R	SALSA MLPA Probemix P065 Marfan Syndrome-1 - 100 reactions
P066-025R	SALSA MLPA Probemix P066 Marfan Syndrome-2 - 25 reactions
P066-050R	SALSA MLPA Probemix P066 Marfan Syndrome-2 - 50 reactions
P066-100R	SALSA MLPA Probemix P066 Marfan Syndrome-2 - 100 reactions
P067-025R	SALSA MLPA Probemix P067 PTCH1 - 25 reactions
P067-050R	SALSA MLPA Probemix P067 PTCH1 - 50 reactions
P067-100R	SALSA MLPA Probemix P067 PTCH1 - 100 reactions
P070-025R	SALSA MLPA Probemix P070 Subtelomeres Mix 2B - 25 reactions
P070-050R	SALSA MLPA Probemix P070 Subtelomeres Mix 2B - 50 reactions
P070-100R	SALSA MLPA Probemix P070 Subtelomeres Mix 2B - 100 reactions
P071-025R	SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3 - 25 reactions
P071-050R	SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3 - 50 reactions
P071-100R	SALSA MLPA Probemix P071 LMNB1-PLP1-NOTCH3 - 100 reactions
P072-025R	SALSA MLPA Probemix P072 MSH6-MUTYH - 25 reactions
P072-050R	SALSA MLPA Probemix P072 MSH6-MUTYH - 50 reactions
P072-100R	SALSA MLPA Probemix P072 MSH6-MUTYH - 100 reactions
P073-025R	SALSA MLPA Probemix P073 IKBKG - 25 reactions
P073-050R	SALSA MLPA Probemix P073 IKBKG - 50 reactions
P073-100R	SALSA MLPA Probemix P073 IKBKG - 100 reactions
P074-025R	SALSA MLPA Probemix P074 AR - 25 reactions
P074-050R	SALSA MLPA Probemix P074 AR - 50 reactions
P074-100R	SALSA MLPA Probemix P074 AR - 100 reactions
P075-025R	SALSA MLPA Probemix P075 TCF4-FOXP1 - 25 reactions
P075-050R	SALSA MLPA Probemix P075 TCF4-FOXP1 - 50 reactions

P075-100R	SALSA MLPA Probemix P075 TCF4-FOXG1 - 100 reactions
P076-025R	SALSA MLPA Probemix P076 ACADVL-SLC22A5 - 25 reactions
P076-050R	SALSA MLPA Probemix P076 ACADVL-SLC22A5 - 50 reactions
P076-100R	SALSA MLPA Probemix P076 ACADVL-SLC22A5 - 100 reactions
P077-025R	SALSA MLPA Probemix P077 BRCA2 Confirmation - 25 reactions
P077-050R	SALSA MLPA Probemix P077 BRCA2 Confirmation - 50 reactions
P077-100R	SALSA MLPA Probemix P077 BRCA2 Confirmation - 100 reactions
P078-025R	SALSA MLPA Probemix P078 Breast tumor - 25 reactions
P078-050R	SALSA MLPA Probemix P078 Breast tumor - 50 reactions
P078-100R	SALSA MLPA Probemix P078 Breast tumor - 100 reactions
P079-025R	SALSA MLPA Probemix P079 OTC - 25 reactions
P079-050R	SALSA MLPA Probemix P079 OTC - 50 reactions
P079-100R	SALSA MLPA Probemix P079 OTC - 100 reactions
P080-025R	SALSA MLPA Probemix P080 Craniofacial - 25 reactions
P080-050R	SALSA MLPA Probemix P080 Craniofacial - 50 reactions
P080-100R	SALSA MLPA Probemix P080 Craniofacial - 100 reactions
P081-025R	SALSA MLPA Probemix P081 NF1 mix 1 - 25 reactions
P081-050R	SALSA MLPA Probemix P081 NF1 mix 1 - 50 reactions
P081-100R	SALSA MLPA Probemix P081 NF1 mix 1 - 100 reactions
P082-025R	SALSA MLPA Probemix P082 NF1 mix 2 - 25 reactions
P082-050R	SALSA MLPA Probemix P082 NF1 mix 2 - 50 reactions
P082-100R	SALSA MLPA Probemix P082 NF1 mix 2 - 100 reactions
P083-025R	SALSA MLPA Probemix P083 CDH1 - 25 reactions
P083-050R	SALSA MLPA Probemix P083 CDH1 - 50 reactions
P083-100R	SALSA MLPA Probemix P083 CDH1 - 100 reactions
P087-025R	SALSA MLPA Probemix P087 BRCA1 Confirmation - 25 reactions
P087-050R	SALSA MLPA Probemix P087 BRCA1 Confirmation - 50 reactions
P087-100R	SALSA MLPA Probemix P087 BRCA1 Confirmation - 100 reactions
P088-025R	SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q - 25 reactions
P088-050R	SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q - 50 reactions
P088-100R	SALSA MLPA Probemix P088 Oligodendroglioma 1p-19q - 100 reactions
P089-025R	SALSA MLPA Probemix P089 TK2 - 25 reactions
P089-050R	SALSA MLPA Probemix P089 TK2 - 50 reactions
P089-100R	SALSA MLPA Probemix P089 TK2 - 100 reactions
P090-025R	SALSA MLPA Probemix P090 BRCA2 - 25 reactions
P090-050R	SALSA MLPA Probemix P090 BRCA2 - 50 reactions
P090-100R	SALSA MLPA Probemix P090 BRCA2 - 100 reactions
P091-025R	SALSA MLPA Probemix P091 CFTR - 25 reactions
P091-050R	SALSA MLPA Probemix P091 CFTR - 50 reactions
P091-100R	SALSA MLPA Probemix P091 CFTR - 100 reactions
P092-025R	SALSA MLPA Probemix P092 ABCC6 - 25 reactions
P092-050R	SALSA MLPA Probemix P092 ABCC6 - 50 reactions
P092-100R	SALSA MLPA Probemix P092 ABCC6 - 100 reactions
P093-025R	SALSA MLPA Probemix P093 HHT/HPAH - 25 reactions
P093-050R	SALSA MLPA Probemix P093 HHT/HPAH - 50 reactions
P093-100R	SALSA MLPA Probemix P093 HHT/HPAH - 100 reactions
P094-025R	SALSA MLPA Probemix P094 MEFV - 25 reactions
P094-050R	SALSA MLPA Probemix P094 MEFV - 50 reactions
P094-100R	SALSA MLPA Probemix P094 MEFV - 100 reactions
P095-025R	SALSA MLPA Probemix P095 Aneuploidy - 25 reactions
P095-050R	SALSA MLPA Probemix P095 Aneuploidy - 50 reactions
P095-100R	SALSA MLPA Probemix P095 Aneuploidy - 100 reactions
P098-025R	SALSA MLPA Probemix P098 Wilson disease - 25 reactions
P098-050R	SALSA MLPA Probemix P098 Wilson disease - 50 reactions

P098-100R	SALSA MLPA Probemix P098 Wilson disease - 100 reactions
P099-025R	SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2 - 25 reactions
P099-050R	SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2 - 50 reactions
P099-100R	SALSA MLPA Probemix P099 GCH1-TH-SGCE-PRRT2 - 100 reactions
P100-025R	SALSA MLPA Probemix P100 MYBPC3 - 25 reactions
P100-050R	SALSA MLPA Probemix P100 MYBPC3 - 50 reactions
P100-100R	SALSA MLPA Probemix P100 MYBPC3 - 100 reactions
P101-025R	SALSA MLPA Probemix P101 STK11 - 25 reactions
P101-050R	SALSA MLPA Probemix P101 STK11 - 50 reactions
P101-100R	SALSA MLPA Probemix P101 STK11 - 100 reactions
P102-025R	SALSA MLPA Probemix P102 HBB - 25 reactions
P102-050R	SALSA MLPA Probemix P102 HBB - 50 reactions
P102-100R	SALSA MLPA Probemix P102 HBB - 100 reactions
P103-025R	SALSA MLPA Probemix P103 DPYD - 25 reactions
P103-050R	SALSA MLPA Probemix P103 DPYD - 50 reactions
P103-100R	SALSA MLPA Probemix P103 DPYD - 100 reactions
P104-025R	SALSA MLPA Probemix P104 Menkes ATP7A - 25 reactions
P104-050R	SALSA MLPA Probemix P104 Menkes ATP7A - 50 reactions
P104-100R	SALSA MLPA Probemix P104 Menkes ATP7A - 100 reactions
P105-025R	SALSA MLPA Probemix P105 Glioma-2 - 25 reactions
P105-050R	SALSA MLPA Probemix P105 Glioma-2 - 50 reactions
P105-100R	SALSA MLPA Probemix P105 Glioma-2 - 100 reactions
P106-025R	SALSA MLPA Probemix P106 X-linked ID - 25 reactions
P106-050R	SALSA MLPA Probemix P106 X-linked ID - 50 reactions
P106-100R	SALSA MLPA Probemix P106 X-linked ID - 100 reactions
P107-025R	SALSA MLPA Probemix P107 Neurometabolic - 25 reactions
P107-050R	SALSA MLPA Probemix P107 Neurometabolic - 50 reactions
P107-100R	SALSA MLPA Probemix P107 Neurometabolic - 100 reactions
P108-025R	SALSA MLPA Probemix P108 SCN5A - 25 reactions
P108-050R	SALSA MLPA Probemix P108 SCN5A - 50 reactions
P108-100R	SALSA MLPA Probemix P108 SCN5A - 100 reactions
P109-025R	SALSA MLPA Probemix P109 ABCB4 - 25 reactions
P109-050R	SALSA MLPA Probemix P109 ABCB4 - 50 reactions
P109-100R	SALSA MLPA Probemix P109 ABCB4 - 100 reactions
P110-025R	SALSA MLPA Probemix P110 FCGR mix 1 - 25 reactions
P110-050R	SALSA MLPA Probemix P110 FCGR mix 1 - 50 reactions
P110-100R	SALSA MLPA Probemix P110 FCGR mix 1 - 100 reactions
P111-025R	SALSA MLPA Probemix P111 FCGR mix 2 - 25 reactions
P111-050R	SALSA MLPA Probemix P111 FCGR mix 2 - 50 reactions
P111-100R	SALSA MLPA Probemix P111 FCGR mix 2 - 100 reactions
P112-025R	SALSA MLPA Probemix P112 PROS1 - 25 reactions
P112-050R	SALSA MLPA Probemix P112 PROS1 - 50 reactions
P112-100R	SALSA MLPA Probemix P112 PROS1 - 100 reactions
P113-025R	SALSA MLPA Probemix P113 FANCB - 25 reactions
P113-050R	SALSA MLPA Probemix P113 FANCB - 50 reactions
P113-100R	SALSA MLPA Probemix P113 FANCB - 100 reactions
P114-025R	SALSA MLPA Probemix P114 Long-QT - 25 reactions
P114-050R	SALSA MLPA Probemix P114 Long-QT - 50 reactions
P114-100R	SALSA MLPA Probemix P114 Long-QT - 100 reactions
P116-025R	SALSA MLPA Probemix P116 SGC - 25 reactions
P116-050R	SALSA MLPA Probemix P116 SGC - 50 reactions
P116-100R	SALSA MLPA Probemix P116 SGC - 100 reactions
P117-025R	SALSA MLPA Probemix P117 ABCC8 - 25 reactions
P117-050R	SALSA MLPA Probemix P117 ABCC8 - 50 reactions

P117-100R	SALSA MLPA Probemix P117 ABCC8 - 100 reactions
P118-025R	SALSA MLPA Probemix P118 WT1 - 25 reactions
P118-050R	SALSA MLPA Probemix P118 WT1 - 50 reactions
P118-100R	SALSA MLPA Probemix P118 WT1 - 100 reactions
P120-025R	SALSA MLPA Probemix P120 PANK2/PLA2G6 - 25 reactions
P120-050R	SALSA MLPA Probemix P120 PANK2/PLA2G6 - 50 reactions
P120-100R	SALSA MLPA Probemix P120 PANK2/PLA2G6 - 100 reactions
P122-025R	SALSA MLPA Probemix P122 NF1-area - 25 reactions
P122-050R	SALSA MLPA Probemix P122 NF1-area - 50 reactions
P122-100R	SALSA MLPA Probemix P122 NF1-area - 100 reactions
P124-025R	SALSA MLPA Probemix P124 TSC1 - 25 reactions
P124-050R	SALSA MLPA Probemix P124 TSC1 - 50 reactions
P124-100R	SALSA MLPA Probemix P124 TSC1 - 100 reactions
P125-025R	SALSA MLPA Probemix P125 Mitochondrial DNA - 25 reactions
P125-050R	SALSA MLPA Probemix P125 Mitochondrial DNA - 50 reactions
P125-100R	SALSA MLPA Probemix P125 Mitochondrial DNA - 100 reactions
P128-025R	SALSA MLPA Probemix P128 CYP450 - 25 reactions
P128-050R	SALSA MLPA Probemix P128 CYP450 - 50 reactions
P128-100R	SALSA MLPA Probemix P128 CYP450 - 100 reactions
P130-025R	SALSA MLPA Probemix P130 CCM mix-1 - 25 reactions
P130-050R	SALSA MLPA Probemix P130 CCM mix-1 - 50 reactions
P130-100R	SALSA MLPA Probemix P130 CCM mix-1 - 100 reactions
P131-025R	SALSA MLPA Probemix P131 CCM mix-2 - 25 reactions
P131-050R	SALSA MLPA Probemix P131 CCM mix-2 - 50 reactions
P131-100R	SALSA MLPA Probemix P131 CCM mix-2 - 100 reactions
P132-025R	SALSA MLPA Probemix P132 Kallmann-1 - 25 reactions
P132-050R	SALSA MLPA Probemix P132 Kallmann-1 - 50 reactions
P132-100R	SALSA MLPA Probemix P132 Kallmann-1 - 100 reactions
P133-025R	SALSA MLPA Probemix P133 Kallmann-2 - 25 reactions
P133-050R	SALSA MLPA Probemix P133 Kallmann-2 - 50 reactions
P133-100R	SALSA MLPA Probemix P133 Kallmann-2 - 100 reactions
P136-025R	SALSA MLPA Probemix P136 Gitelman Syndrome - 25 reactions
P136-050R	SALSA MLPA Probemix P136 Gitelman Syndrome - 50 reactions
P136-100R	SALSA MLPA Probemix P136 Gitelman Syndrome - 100 reactions
P137-025R	SALSA MLPA Probemix P137 SCN1A - 25 reactions
P137-050R	SALSA MLPA Probemix P137 SCN1A - 50 reactions
P137-100R	SALSA MLPA Probemix P137 SCN1A - 100 reactions
P138-025R	SALSA MLPA Probemix P138 SLC2A1-STXBP1 - 25 reactions
P138-050R	SALSA MLPA Probemix P138 SLC2A1-STXBP1 - 50 reactions
P138-100R	SALSA MLPA Probemix P138 SLC2A1-STXBP1 - 100 reactions
P140-025R	SALSA MLPA Probemix P140 HBA - 25 reactions
P140-050R	SALSA MLPA Probemix P140 HBA - 50 reactions
P140-100R	SALSA MLPA Probemix P140 HBA - 100 reactions
P141-025R	SALSA MLPA Probemix P141 NIPBL-1 - 25 reactions
P141-050R	SALSA MLPA Probemix P141 NIPBL-1 - 50 reactions
P141-100R	SALSA MLPA Probemix P141 NIPBL-1 - 100 reactions
P142-025R	SALSA MLPA Probemix P142 NIPBL-2 - 25 reactions
P142-050R	SALSA MLPA Probemix P142 NIPBL-2 - 50 reactions
P142-100R	SALSA MLPA Probemix P142 NIPBL-2 - 100 reactions
P143-025R	SALSA MLPA Probemix P143 MFN2-MPZ - 25 reactions
P143-050R	SALSA MLPA Probemix P143 MFN2-MPZ - 50 reactions
P143-100R	SALSA MLPA Probemix P143 MFN2-MPZ - 100 reactions
P147-025R	SALSA MLPA Probemix P147 1p36 - 25 reactions
P147-050R	SALSA MLPA Probemix P147 1p36 - 50 reactions

P147-100R	SALSA MLPA Probemix P147 1p36 - 100 reactions
P148-025R	SALSA MLPA Probemix P148 TGFBR1-TGFBR2 - 25 reactions
P148-050R	SALSA MLPA Probemix P148 TGFBR1-TGFBR2 - 50 reactions
P148-100R	SALSA MLPA Probemix P148 TGFBR1-TGFBR2 - 100 reactions
P151-025R	SALSA MLPA Probemix P151 ABCA4 mix-1 - 25 reactions
P151-050R	SALSA MLPA Probemix P151 ABCA4 mix-1 - 50 reactions
P151-100R	SALSA MLPA Probemix P151 ABCA4 mix-1 - 100 reactions
P152-025R	SALSA MLPA Probemix P152 ABCA4 mix-2 - 25 reactions
P152-050R	SALSA MLPA Probemix P152 ABCA4 mix-2 - 50 reactions
P152-100R	SALSA MLPA Probemix P152 ABCA4 mix-2 - 100 reactions
P153-025R	SALSA MLPA Probemix P153 EYA1 - 25 reactions
P153-050R	SALSA MLPA Probemix P153 EYA1 - 50 reactions
P153-100R	SALSA MLPA Probemix P153 EYA1 - 100 reactions
P154-025R	SALSA MLPA Probemix P154 GPC3-GPC4 - 25 reactions
P154-050R	SALSA MLPA Probemix P154 GPC3-GPC4 - 50 reactions
P154-100R	SALSA MLPA Probemix P154 GPC3-GPC4 - 100 reactions
P155-025R	SALSA MLPA Probemix P155 EDS - 25 reactions
P155-050R	SALSA MLPA Probemix P155 EDS - 50 reactions
P155-100R	SALSA MLPA Probemix P155 EDS - 100 reactions
P156-025R	SALSA MLPA Probemix P156 GALT - 25 reactions
P156-050R	SALSA MLPA Probemix P156 GALT - 50 reactions
P156-100R	SALSA MLPA Probemix P156 GALT - 100 reactions
P158-025R	SALSA MLPA Probemix P158 JPS - 25 reactions
P158-050R	SALSA MLPA Probemix P158 JPS - 50 reactions
P158-100R	SALSA MLPA Probemix P158 JPS - 100 reactions
P159-025R	SALSA MLPA Probemix P159 GLA - 25 reactions
P159-050R	SALSA MLPA Probemix P159 GLA - 50 reactions
P159-100R	SALSA MLPA Probemix P159 GLA - 100 reactions
P160-025R	SALSA MLPA Probemix P160 STS - 25 reactions
P160-050R	SALSA MLPA Probemix P160 STS - 50 reactions
P160-100R	SALSA MLPA Probemix P160 STS - 100 reactions
P163-025R	SALSA MLPA Probemix P163 GJB-WFS1-POU3F4 - 25 reactions
P163-050R	SALSA MLPA Probemix P163 GJB-WFS1-POU3F4 - 50 reactions
P163-100R	SALSA MLPA Probemix P163 GJB-WFS1-POU3F4 - 100 reactions
P164-025R	SALSA MLPA Probemix P164 IDS - 25 reactions
P164-050R	SALSA MLPA Probemix P164 IDS - 50 reactions
P164-100R	SALSA MLPA Probemix P164 IDS - 100 reactions
P165-025R	SALSA MLPA Probemix P165 HSP mix-1 - 25 reactions
P165-050R	SALSA MLPA Probemix P165 HSP mix-1 - 50 reactions
P165-100R	SALSA MLPA Probemix P165 HSP mix-1 - 100 reactions
P166-025R	SALSA MLPA Probemix P166 KCNQ2 - 25 reactions
P166-050R	SALSA MLPA Probemix P166 KCNQ2 - 50 reactions
P166-100R	SALSA MLPA Probemix P166 KCNQ2 - 100 reactions
P168-025R	SALSA MLPA Probemix P168 ARVC-PKP2 - 25 reactions
P168-050R	SALSA MLPA Probemix P168 ARVC-PKP2 - 50 reactions
P168-100R	SALSA MLPA Probemix P168 ARVC-PKP2 - 100 reactions
P169-025R	SALSA MLPA Probemix P169 Hirschsprung-1 - 25 reactions
P169-050R	SALSA MLPA Probemix P169 Hirschsprung-1 - 50 reactions
P169-100R	SALSA MLPA Probemix P169 Hirschsprung-1 - 100 reactions
P170-025R	SALSA MLPA Probemix P170 APP - 25 reactions
P170-050R	SALSA MLPA Probemix P170 APP - 50 reactions
P170-100R	SALSA MLPA Probemix P170 APP - 100 reactions
P175-025R	SALSA MLPA Probemix P175 Tumour Gain - 25 reactions
P175-050R	SALSA MLPA Probemix P175 Tumour Gain - 50 reactions

P175-100R	SALSA MLPA Probemix P175 Tumour Gain - 100 reactions
P176-025R	SALSA MLPA Probemix P176 CAPN3 - 25 reactions
P176-050R	SALSA MLPA Probemix P176 CAPN3 - 50 reactions
P176-100R	SALSA MLPA Probemix P176 CAPN3 - 100 reactions
P177-025R	SALSA MLPA Probemix P177 CASR - 25 reactions
P177-050R	SALSA MLPA Probemix P177 CASR - 50 reactions
P177-100R	SALSA MLPA Probemix P177 CASR - 100 reactions
P178-025R	SALSA MLPA Probemix P178 F8 - 25 reactions
P178-050R	SALSA MLPA Probemix P178 F8 - 50 reactions
P178-100R	SALSA MLPA Probemix P178 F8 - 100 reactions
P179-025R	SALSA MLPA Probemix P179 Limb-1 - 25 reactions
P179-050R	SALSA MLPA Probemix P179 Limb-1 - 50 reactions
P179-100R	SALSA MLPA Probemix P179 Limb-1 - 100 reactions
P180-025R	SALSA MLPA Probemix P180 Limb-2 - 25 reactions
P180-050R	SALSA MLPA Probemix P180 Limb-2 - 50 reactions
P180-100R	SALSA MLPA Probemix P180 Limb-2 - 100 reactions
P181-025R	SALSA MLPA Probemix P181 Centromere mix 1 - 25 reactions
P181-050R	SALSA MLPA Probemix P181 Centromere mix 1 - 50 reactions
P181-100R	SALSA MLPA Probemix P181 Centromere mix 1 - 100 reactions
P182-025R	SALSA MLPA Probemix P182 Centromere mix 2 - 25 reactions
P182-050R	SALSA MLPA Probemix P182 Centromere mix 2 - 50 reactions
P182-100R	SALSA MLPA Probemix P182 Centromere mix 2 - 100 reactions
P183-025R	SALSA MLPA Probemix P183 EDA-EDAR-EDARADD - 25 reactions
P183-050R	SALSA MLPA Probemix P183 EDA-EDAR-EDARADD - 50 reactions
P183-100R	SALSA MLPA Probemix P183 EDA-EDAR-EDARADD - 100 reactions
P184-025R	SALSA MLPA Probemix P184 JAG1 - 25 reactions
P184-050R	SALSA MLPA Probemix P184 JAG1 - 50 reactions
P184-100R	SALSA MLPA Probemix P184 JAG1 - 100 reactions
P185-025R	SALSA MLPA Probemix P185 Intersex - 25 reactions
P185-050R	SALSA MLPA Probemix P185 Intersex - 50 reactions
P185-100R	SALSA MLPA Probemix P185 Intersex - 100 reactions
P186-025R	SALSA MLPA Probemix P186 PAX3-MITF-SOX10 - 25 reactions
P186-050R	SALSA MLPA Probemix P186 PAX3-MITF-SOX10 - 50 reactions
P186-100R	SALSA MLPA Probemix P186 PAX3-MITF-SOX10 - 100 reactions
P187-025R	SALSA MLPA Probemix P187 HPE - 25 reactions
P187-050R	SALSA MLPA Probemix P187 HPE - 50 reactions
P187-100R	SALSA MLPA Probemix P187 HPE - 100 reactions
P188-025R	SALSA MLPA Probemix P188 22q13 - 25 reactions
P188-050R	SALSA MLPA Probemix P188 22q13 - 50 reactions
P188-100R	SALSA MLPA Probemix P188 22q13 - 100 reactions
P189-025R	SALSA MLPA Probemix P189 CDKL5/ARX/FOXP1 - 25 reactions
P189-050R	SALSA MLPA Probemix P189 CDKL5/ARX/FOXP1 - 50 reactions
P189-100R	SALSA MLPA Probemix P189 CDKL5/ARX/FOXP1 - 100 reactions
P190-025R	SALSA MLPA Probemix P190 CHEK2 - 25 reactions
P190-050R	SALSA MLPA Probemix P190 CHEK2 - 50 reactions
P190-100R	SALSA MLPA Probemix P190 CHEK2 - 100 reactions
P191-025R	SALSA MLPA Probemix P191 Alport-mix 1 - 25 reactions
P191-050R	SALSA MLPA Probemix P191 Alport-mix 1 - 50 reactions
P191-100R	SALSA MLPA Probemix P191 Alport-mix 1 - 100 reactions
P192-025R	SALSA MLPA Probemix P192 Alport-mix 2 - 25 reactions
P192-050R	SALSA MLPA Probemix P192 Alport-mix 2 - 50 reactions
P192-100R	SALSA MLPA Probemix P192 Alport-mix 2 - 100 reactions
P193-025R	SALSA MLPA Probemix P193 NPC1-NPC2-SMPD1 - 25 reactions
P193-050R	SALSA MLPA Probemix P193 NPC1-NPC2-SMPD1 - 50 reactions

P193-100R	SALSA MLPA Probemix P193 NPC1-NPC2-SMPD1 - 100 reactions
P196-025R	SALSA MLPA Probemix P196 TNNT2-BAG3 - 25 reactions
P196-050R	SALSA MLPA Probemix P196 TNNT2-BAG3 - 50 reactions
P196-100R	SALSA MLPA Probemix P196 TNNT2-BAG3 - 100 reactions
P197-025R	SALSA MLPA Probemix P197 KCNQ3 - 25 reactions
P197-050R	SALSA MLPA Probemix P197 KCNQ3 - 50 reactions
P197-100R	SALSA MLPA Probemix P197 KCNQ3 - 100 reactions
P198-025R	SALSA MLPA Probemix P198 FH - 25 reactions
P198-050R	SALSA MLPA Probemix P198 FH - 50 reactions
P198-100R	SALSA MLPA Probemix P198 FH - 100 reactions
P199-025R	SALSA MLPA Probemix P199 HEXA - 25 reactions
P199-050R	SALSA MLPA Probemix P199 HEXA - 50 reactions
P199-100R	SALSA MLPA Probemix P199 HEXA - 100 reactions
P200-100R	SALSA MLPA Probemix P200 Reference-1 - 100 reactions
P201-025R	SALSA MLPA Probemix P201 CHARGE - 25 reactions
P201-050R	SALSA MLPA Probemix P201 CHARGE - 50 reactions
P201-100R	SALSA MLPA Probemix P201 CHARGE - 100 reactions
P202-025R	SALSA MLPA Probemix P202 IKZF1-ERG - 25 reactions
P202-050R	SALSA MLPA Probemix P202 IKZF1-ERG - 50 reactions
P202-100R	SALSA MLPA Probemix P202 IKZF1-ERG - 100 reactions
P203-025R	SALSA MLPA Probemix P203 PKLR - 25 reactions
P203-050R	SALSA MLPA Probemix P203 PKLR - 50 reactions
P203-100R	SALSA MLPA Probemix P203 PKLR - 100 reactions
P205-025R	SALSA MLPA Probemix P205 SH2D1A-XIAP-ITK - 25 reactions
P205-050R	SALSA MLPA Probemix P205 SH2D1A-XIAP-ITK - 50 reactions
P205-100R	SALSA MLPA Probemix P205 SH2D1A-XIAP-ITK - 100 reactions
P207-025R	SALSA MLPA Probemix P207 F9 - 25 reactions
P207-050R	SALSA MLPA Probemix P207 F9 - 50 reactions
P207-100R	SALSA MLPA Probemix P207 F9 - 100 reactions
P209-025R	SALSA MLPA Probemix P209 Glycine Encephalopathy - 25 reactions
P209-050R	SALSA MLPA Probemix P209 Glycine Encephalopathy - 50 reactions
P209-100R	SALSA MLPA Probemix P209 Glycine Encephalopathy - 100 reactions
P210-025R	SALSA MLPA Probemix P210 BTK - 25 reactions
P210-050R	SALSA MLPA Probemix P210 BTK - 50 reactions
P210-100R	SALSA MLPA Probemix P210 BTK - 100 reactions
P211-025R	SALSA MLPA Probemix P211 HSP region - 25 reactions
P211-050R	SALSA MLPA Probemix P211 HSP region - 50 reactions
P211-100R	SALSA MLPA Probemix P211 HSP region - 100 reactions
P212-025R	SALSA MLPA Probemix P212 DBA - 25 reactions
P212-050R	SALSA MLPA Probemix P212 DBA - 50 reactions
P212-100R	SALSA MLPA Probemix P212 DBA - 100 reactions
P213-025R	SALSA MLPA Probemix P213 HSP mix-2 - 25 reactions
P213-050R	SALSA MLPA Probemix P213 HSP mix-2 - 50 reactions
P213-100R	SALSA MLPA Probemix P213 HSP mix-2 - 100 reactions
P214-025R	SALSA MLPA Probemix P214 COL2A1 - 25 reactions
P214-050R	SALSA MLPA Probemix P214 COL2A1 - 50 reactions
P214-100R	SALSA MLPA Probemix P214 COL2A1 - 100 reactions
P215-025R	SALSA MLPA Probemix P215 EXT - 25 reactions
P215-050R	SALSA MLPA Probemix P215 EXT - 50 reactions
P215-100R	SALSA MLPA Probemix P215 EXT - 100 reactions
P216-025R	SALSA MLPA Probemix P216 GHD mix1 - 25 reactions
P216-050R	SALSA MLPA Probemix P216 GHD mix1 - 50 reactions
P216-100R	SALSA MLPA Probemix P216 GHD mix1 - 100 reactions
P217-025R	SALSA MLPA Probemix P217 IGF1R - 25 reactions

P217-050R	SALSA MLPA Probemix P217 IGF1R - 50 reactions
P217-100R	SALSA MLPA Probemix P217 IGF1R - 100 reactions
P218-025R	SALSA MLPA Probemix P218 LPL - 25 reactions
P218-050R	SALSA MLPA Probemix P218 LPL - 50 reactions
P218-100R	SALSA MLPA Probemix P218 LPL - 100 reactions
P219-025R	SALSA MLPA Probemix P219 PAX6 - 25 reactions
P219-050R	SALSA MLPA Probemix P219 PAX6 - 50 reactions
P219-100R	SALSA MLPA Probemix P219 PAX6 - 100 reactions
P220-025R	SALSA MLPA Probemix P220 Obesity - 25 reactions
P220-050R	SALSA MLPA Probemix P220 Obesity - 50 reactions
P220-100R	SALSA MLPA Probemix P220 Obesity - 100 reactions
P221-025R	SALSA MLPA Probemix P221 LCA mix-1 - 25 reactions
P221-050R	SALSA MLPA Probemix P221 LCA mix-1 - 50 reactions
P221-100R	SALSA MLPA Probemix P221 LCA mix-1 - 100 reactions
P222-025R	SALSA MLPA Probemix P222 LCA mix-2 - 25 reactions
P222-050R	SALSA MLPA Probemix P222 LCA mix-2 - 50 reactions
P222-100R	SALSA MLPA Probemix P222 LCA mix-2 - 100 reactions
P223-025R	SALSA MLPA Probemix P223 PHEX - 25 reactions
P223-050R	SALSA MLPA Probemix P223 PHEX - 50 reactions
P223-100R	SALSA MLPA Probemix P223 PHEX - 100 reactions
P225-025R	SALSA MLPA Probemix P225 PTEN - 25 reactions
P225-050R	SALSA MLPA Probemix P225 PTEN - 50 reactions
P225-100R	SALSA MLPA Probemix P225 PTEN - 100 reactions
P226-025R	SALSA MLPA Probemix P226 SDH - 25 reactions
P226-050R	SALSA MLPA Probemix P226 SDH - 50 reactions
P226-100R	SALSA MLPA Probemix P226 SDH - 100 reactions
P227-025R	SALSA MLPA Probemix P227 SerpinC1 - 25 reactions
P227-050R	SALSA MLPA Probemix P227 SerpinC1 - 50 reactions
P227-100R	SALSA MLPA Probemix P227 SerpinC1 - 100 reactions
P229-025R	SALSA MLPA Probemix P229 OPA1 - 25 reactions
P229-050R	SALSA MLPA Probemix P229 OPA1 - 50 reactions
P229-100R	SALSA MLPA Probemix P229 OPA1 - 100 reactions
P232-025R	SALSA MLPA Probemix P232 FGD1 - 25 reactions
P232-050R	SALSA MLPA Probemix P232 FGD1 - 50 reactions
P232-100R	SALSA MLPA Probemix P232 FGD1 - 100 reactions
P234-025R	SALSA MLPA Probemix P234 GATA3 - GATA4 - 25 reactions
P234-050R	SALSA MLPA Probemix P234 GATA3 - GATA4 - 50 reactions
P234-100R	SALSA MLPA Probemix P234 GATA3 - GATA4 - 100 reactions
P235-025R	SALSA MLPA Probemix P235 Retinitis - 25 reactions
P235-050R	SALSA MLPA Probemix P235 Retinitis - 50 reactions
P235-100R	SALSA MLPA Probemix P235 Retinitis - 100 reactions
P236-025R	SALSA MLPA Probemix P236 CFH Region - 25 reactions
P236-050R	SALSA MLPA Probemix P236 CFH Region - 50 reactions
P236-100R	SALSA MLPA Probemix P236 CFH Region - 100 reactions
P237-025R	SALSA MLPA Probemix P237 DNAI1 - 25 reactions
P237-050R	SALSA MLPA Probemix P237 DNAI1 - 50 reactions
P237-100R	SALSA MLPA Probemix P237 DNAI1 - 100 reactions
P238-025R	SALSA MLPA Probemix P238 DNAH5 - 25 reactions
P238-050R	SALSA MLPA Probemix P238 DNAH5 - 50 reactions
P238-100R	SALSA MLPA Probemix P238 DNAH5 - 100 reactions
P239-025R	SALSA MLPA Probemix P239 BRCA1 region - 25 reactions
P239-050R	SALSA MLPA Probemix P239 BRCA1 region - 50 reactions
P239-100R	SALSA MLPA Probemix P239 BRCA1 region - 100 reactions
P240-025R	SALSA MLPA Probemix P240 BRIP1/CHEK1 - 25 reactions

P240-050R	SALSA MLPA Probemix P240 BRIP1/CHEK1 - 50 reactions
P240-100R	SALSA MLPA Probemix P240 BRIP1/CHEK1 - 100 reactions
P241-025R	SALSA MLPA Probemix P241 MODY Mix 1 - 25 reactions
P241-050R	SALSA MLPA Probemix P241 MODY Mix 1 - 50 reactions
P241-100R	SALSA MLPA Probemix P241 MODY Mix 1 - 100 reactions
P242-025R	SALSA MLPA Probemix P242 Pancreatitis - 25 reactions
P242-050R	SALSA MLPA Probemix P242 Pancreatitis - 50 reactions
P242-100R	SALSA MLPA Probemix P242 Pancreatitis - 100 reactions
P243-025R	SALSA MLPA Probemix P243 SERPING1 - F12 - 25 reactions
P243-050R	SALSA MLPA Probemix P243 SERPING1 - F12 - 50 reactions
P243-100R	SALSA MLPA Probemix P243 SERPING1 - F12 - 100 reactions
P244-025R	SALSA MLPA Probemix P244 AIP-MEN1-CDKN1B - 25 reactions
P244-050R	SALSA MLPA Probemix P244 AIP-MEN1-CDKN1B - 50 reactions
P244-100R	SALSA MLPA Probemix P244 AIP-MEN1-CDKN1B - 100 reactions
P245-025R	SALSA MLPA Probemix P245 Microdeletion Syndromes-1A - 25 reactions
P245-050R	SALSA MLPA Probemix P245 Microdeletion Syndromes-1A - 50 reactions
P245-100R	SALSA MLPA Probemix P245 Microdeletion Syndromes-1A - 100 reactions
P248-025R	SALSA MLPA Probemix P248 MLH1-MSH2 Confirmation - 25 reactions
P248-050R	SALSA MLPA Probemix P248 MLH1-MSH2 Confirmation - 50 reactions
P248-100R	SALSA MLPA Probemix P248 MLH1-MSH2 Confirmation - 100 reactions
P250-025R	SALSA MLPA Probemix P250 DiGeorge - 25 reactions
P250-050R	SALSA MLPA Probemix P250 DiGeorge - 50 reactions
P250-100R	SALSA MLPA Probemix P250 DiGeorge - 100 reactions
P251-025R	SALSA MLPA Probemix P251 NB mix 1 - 25 reactions
P251-050R	SALSA MLPA Probemix P251 NB mix 1 - 50 reactions
P252-025R	SALSA MLPA Probemix P252 NB mix 2 - 25 reactions
P252-050R	SALSA MLPA Probemix P252 NB mix 2 - 50 reactions
P252-100R	SALSA MLPA Probemix P252 NB mix 2 - 100 reactions
P253-025R	SALSA MLPA Probemix P253 NB mix 3 - 25 reactions
P253-050R	SALSA MLPA Probemix P253 NB mix 3 - 50 reactions
P254-025R	SALSA MLPA Probemix P254 PSEN1 - 25 reactions
P254-050R	SALSA MLPA Probemix P254 PSEN1 - 50 reactions
P254-100R	SALSA MLPA Probemix P254 PSEN1 - 100 reactions
P255-025R	SALSA MLPA Probemix P255 ALDOB-FBP1 - 25 reactions
P255-050R	SALSA MLPA Probemix P255 ALDOB-FBP1 - 50 reactions
P255-100R	SALSA MLPA Probemix P255 ALDOB-FBP1 - 100 reactions
P256-025R	SALSA MLPA Probemix P256 FLCN - 25 reactions
P256-050R	SALSA MLPA Probemix P256 FLCN - 50 reactions
P256-100R	SALSA MLPA Probemix P256 FLCN - 100 reactions
P257-025R	SALSA MLPA Probemix P257 TERT-DKC1 - 25 reactions
P257-050R	SALSA MLPA Probemix P257 TERT-DKC1 - 50 reactions
P257-100R	SALSA MLPA Probemix P257 TERT-DKC1 - 100 reactions
P258-025R	SALSA MLPA Probemix P258 SMARCB1 - 25 reactions
P258-050R	SALSA MLPA Probemix P258 SMARCB1 - 50 reactions
P258-100R	SALSA MLPA Probemix P258 SMARCB1 - 100 reactions
P260-025R	SALSA MLPA Probemix P260 PALB2-RAD50-RAD51C-RAD51D - 25 reactions
P260-050R	SALSA MLPA Probemix P260 PALB2-RAD50-RAD51C-RAD51D - 50 reactions
P260-100R	SALSA MLPA Probemix P260 PALB2-RAD50-RAD51C-RAD51D - 100 reactions
P262-025R	SALSA MLPA Probemix P262 GHI - 25 reactions
P262-050R	SALSA MLPA Probemix P262 GHI - 50 reactions
P262-100R	SALSA MLPA Probemix P262 GHI - 100 reactions
P265-025R	SALSA MLPA Probemix P265 PROC - 25 reactions
P265-050R	SALSA MLPA Probemix P265 PROC - 50 reactions

P265-100R	SALSA MLPA Probemix P265 PROC - 100 reactions
P266-025R	SALSA MLPA Probemix P266 CLCNKB - 25 reactions
P266-050R	SALSA MLPA Probemix P266 CLCNKB - 50 reactions
P266-100R	SALSA MLPA Probemix P266 CLCNKB - 100 reactions
P268-025R	SALSA MLPA Probemix P268 DYSF - 25 reactions
P268-050R	SALSA MLPA Probemix P268 DYSF - 50 reactions
P268-100R	SALSA MLPA Probemix P268 DYSF - 100 reactions
P271-025R	SALSA MLPA Probemix P271 COL1A1 - 25 reactions
P271-050R	SALSA MLPA Probemix P271 COL1A1 - 50 reactions
P271-100R	SALSA MLPA Probemix P271 COL1A1 - 100 reactions
P272-025R	SALSA MLPA Probemix P272 COL1A2 - 25 reactions
P272-050R	SALSA MLPA Probemix P272 COL1A2 - 50 reactions
P272-100R	SALSA MLPA Probemix P272 COL1A2 - 100 reactions
P275-025R	SALSA MLPA Probemix P275 MAPT-GRN - 25 reactions
P275-050R	SALSA MLPA Probemix P275 MAPT-GRN - 50 reactions
P275-100R	SALSA MLPA Probemix P275 MAPT-GRN - 100 reactions
P278-025R	SALSA MLPA Probemix P278 PCCA-PCCB - 25 reactions
P278-050R	SALSA MLPA Probemix P278 PCCA-PCCB - 50 reactions
P278-100R	SALSA MLPA Probemix P278 PCCA-PCCB - 100 reactions
P279-025R	SALSA MLPA Probemix P279 CACNA1A - 25 reactions
P279-050R	SALSA MLPA Probemix P279 CACNA1A - 50 reactions
P279-100R	SALSA MLPA Probemix P279 CACNA1A - 100 reactions
P280-025R	SALSA MLPA Probemix P280 SLC26A4 - 25 reactions
P280-050R	SALSA MLPA Probemix P280 SLC26A4 - 50 reactions
P280-100R	SALSA MLPA Probemix P280 SLC26A4 - 100 reactions
P283-025R	SALSA MLPA Probemix P283 TPMT - 25 reactions
P283-050R	SALSA MLPA Probemix P283 TPMT - 50 reactions
P283-100R	SALSA MLPA Probemix P283 TPMT - 100 reactions
P285-025R	SALSA MLPA Probemix P285 LRP5 - 25 reactions
P285-050R	SALSA MLPA Probemix P285 LRP5 - 50 reactions
P285-100R	SALSA MLPA Probemix P285 LRP5 - 100 reactions
P289-025R	SALSA MLPA Probemix P289 LMX1B - 25 reactions
P289-050R	SALSA MLPA Probemix P289 LMX1B - 50 reactions
P289-100R	SALSA MLPA Probemix P289 LMX1B - 100 reactions
P292-025R	SALSA MLPA Probemix P292 PCDH15 - 25 reactions
P292-050R	SALSA MLPA Probemix P292 PCDH15 - 50 reactions
P292-100R	SALSA MLPA Probemix P292 PCDH15 - 100 reactions
P294-025R	SALSA MLPA Probemix P294 Tumour Loss - 25 reactions
P294-050R	SALSA MLPA Probemix P294 Tumour Loss - 50 reactions
P294-100R	SALSA MLPA Probemix P294 Tumour Loss - 100 reactions
P295-025R	SALSA MLPA Probemix P295 SPRED1 - 25 reactions
P295-050R	SALSA MLPA Probemix P295 SPRED1 - 50 reactions
P295-100R	SALSA MLPA Probemix P295 SPRED1 - 100 reactions
P296-025R	SALSA MLPA Probemix P296 aHUS - 25 reactions
P296-050R	SALSA MLPA Probemix P296 aHUS - 50 reactions
P296-100R	SALSA MLPA Probemix P296 aHUS - 100 reactions
P297-025R	SALSA MLPA Probemix P297 Microdeletion-2 - 25 reactions
P297-050R	SALSA MLPA Probemix P297 Microdeletion-2 - 50 reactions
P297-100R	SALSA MLPA Probemix P297 Microdeletion-2 - 100 reactions
P298-025R	SALSA MLPA Probemix P298 BRAF-HRAS-KRAS-NRAS - 25 reactions
P298-050R	SALSA MLPA Probemix P298 BRAF-HRAS-KRAS-NRAS - 50 reactions
P298-100R	SALSA MLPA Probemix P298 BRAF-HRAS-KRAS-NRAS - 100 reactions
P300-100R	SALSA MLPA Probemix P300 Reference-2 - 100 reactions
P301-025R	SALSA MLPA Probemix P301 Medulloblastoma mix 1 - 25 reactions

P301-050R	SALSA MLPA Probemix P301 Medulloblastoma mix 1 - 50 reactions
P302-025R	SALSA MLPA Probemix P302 Medulloblastoma mix 2 - 25 reactions
P302-050R	SALSA MLPA Probemix P302 Medulloblastoma mix 2 - 50 reactions
P303-025R	SALSA MLPA Probemix P303 Medulloblastoma mix 3 - 25 reactions
P303-050R	SALSA MLPA Probemix P303 Medulloblastoma mix 3 - 50 reactions
P305-025R	SALSA MLPA Probemix P305 AGXT - 25 reactions
P305-050R	SALSA MLPA Probemix P305 AGXT - 50 reactions
P305-100R	SALSA MLPA Probemix P305 AGXT - 100 reactions
P306-025R	SALSA MLPA Probemix P306 SPG11 - 25 reactions
P306-050R	SALSA MLPA Probemix P306 SPG11 - 50 reactions
P306-100R	SALSA MLPA Probemix P306 SPG11 - 100 reactions
P307-025R	SALSA MLPA Probemix P307 SEPT9 - 25 reactions
P307-050R	SALSA MLPA Probemix P307 SEPT9 - 50 reactions
P307-100R	SALSA MLPA Probemix P307 SEPT9 - 100 reactions
P308-025R	SALSA MLPA Probemix P308 MET - 25 reactions
P308-050R	SALSA MLPA Probemix P308 MET - 50 reactions
P308-100R	SALSA MLPA Probemix P308 MET - 100 reactions
P309-025R	SALSA MLPA Probemix P309 MTM1 - 25 reactions
P309-050R	SALSA MLPA Probemix P309 MTM1 - 50 reactions
P309-100R	SALSA MLPA Probemix P309 MTM1 - 100 reactions
P310-025R	SALSA MLPA Probemix P310 TCOF1 - 25 reactions
P310-050R	SALSA MLPA Probemix P310 TCOF1 - 50 reactions
P310-100R	SALSA MLPA Probemix P310 TCOF1 - 100 reactions
P311-025R	SALSA MLPA Probemix P311 CHD - 25 reactions
P311-050R	SALSA MLPA Probemix P311 CHD - 50 reactions
P311-100R	SALSA MLPA Probemix P311 CHD - 100 reactions
P312-025R	SALSA MLPA Probemix P312 POR - 25 reactions
P312-050R	SALSA MLPA Probemix P312 POR - 50 reactions
P312-100R	SALSA MLPA Probemix P312 POR - 100 reactions
P313-025R	SALSA MLPA Probemix P313 CREBBP - 25 reactions
P313-050R	SALSA MLPA Probemix P313 CREBBP - 50 reactions
P313-100R	SALSA MLPA Probemix P313 CREBBP - 100 reactions
P314-025R	SALSA MLPA Probemix P314 ABCA3-SFTPC - 25 reactions
P314-050R	SALSA MLPA Probemix P314 ABCA3-SFTPC - 50 reactions
P314-100R	SALSA MLPA Probemix P314 ABCA3-SFTPC - 100 reactions
P315-025R	SALSA MLPA Probemix P315 EGFR - 25 reactions
P315-050R	SALSA MLPA Probemix P315 EGFR - 50 reactions
P315-100R	SALSA MLPA Probemix P315 EGFR - 100 reactions
P316-025R	SALSA MLPA Probemix P316 Recessive Ataxias - 25 reactions
P316-050R	SALSA MLPA Probemix P316 Recessive Ataxias - 50 reactions
P316-100R	SALSA MLPA Probemix P316 Recessive Ataxias - 100 reactions
P318-025R	SALSA MLPA Probemix P318 Hirschsprung-2 - 25 reactions
P318-050R	SALSA MLPA Probemix P318 Hirschsprung-2 - 50 reactions
P318-100R	SALSA MLPA Probemix P318 Hirschsprung-2 - 100 reactions
P319-025R	SALSA MLPA Probemix P319 Thyroid - 25 reactions
P319-050R	SALSA MLPA Probemix P319 Thyroid - 50 reactions
P319-100R	SALSA MLPA Probemix P319 Thyroid - 100 reactions
P321-025R	SALSA MLPA Probemix P321 VPS13B mix 1 - 25 reactions
P321-050R	SALSA MLPA Probemix P321 VPS13B mix 1 - 50 reactions
P321-100R	SALSA MLPA Probemix P321 VPS13B mix 1 - 100 reactions
P322-025R	SALSA MLPA Probemix P322 VPS13B mix 2 - 25 reactions
P322-050R	SALSA MLPA Probemix P322 VPS13B mix 2 - 50 reactions
P322-100R	SALSA MLPA Probemix P322 VPS13B mix 2 - 100 reactions
P323-025R	SALSA MLPA Probemix P323 CDK4-HMGA2-MDM2 - 25 reactions

P323-050R	SALSA MLPA Probemix P323 CDK4-HMGA2-MDM2 - 50 reactions
P323-100R	SALSA MLPA Probemix P323 CDK4-HMGA2-MDM2 - 100 reactions
P324-025R	SALSA MLPA Probemix P324 22q11 - 25 reactions
P324-050R	SALSA MLPA Probemix P324 22q11 - 50 reactions
P324-100R	SALSA MLPA Probemix P324 22q11 - 100 reactions
P325-025R	SALSA MLPA Probemix P325 OCA2 - 25 reactions
P325-050R	SALSA MLPA Probemix P325 OCA2 - 50 reactions
P325-100R	SALSA MLPA Probemix P325 OCA2 - 100 reactions
P326-025R	SALSA MLPA Probemix P326 LARGE1 - 25 reactions
P326-050R	SALSA MLPA Probemix P326 LARGE1 - 50 reactions
P326-100R	SALSA MLPA Probemix P326 LARGE1 - 100 reactions
P327-025R	SALSA MLPA Probemix P327 iAMP21-ERG - 25 reactions
P327-050R	SALSA MLPA Probemix P327 iAMP21-ERG - 50 reactions
P327-100R	SALSA MLPA Probemix P327 iAMP21-ERG - 100 reactions
P328-025R	SALSA MLPA Probemix P328 EYS - 25 reactions
P328-050R	SALSA MLPA Probemix P328 EYS - 50 reactions
P328-100R	SALSA MLPA Probemix P328 EYS - 100 reactions
P329-025R	SALSA MLPA Probemix P329 CRLF2-CSF2RA-IL3RA - 25 reactions
P329-050R	SALSA MLPA Probemix P329 CRLF2-CSF2RA-IL3RA - 50 reactions
P329-100R	SALSA MLPA Probemix P329 CRLF2-CSF2RA-IL3RA - 100 reactions
P330-025R	SALSA MLPA Probemix P330 PCDH19 - 25 reactions
P330-050R	SALSA MLPA Probemix P330 PCDH19 - 50 reactions
P330-100R	SALSA MLPA Probemix P330 PCDH19 - 100 reactions
P331-025R	SALSA MLPA Probemix P331 COL5A1 MIX-1 - 25 reactions
P331-050R	SALSA MLPA Probemix P331 COL5A1 MIX-1 - 50 reactions
P331-100R	SALSA MLPA Probemix P331 COL5A1 MIX-1 - 100 reactions
P332-025R	SALSA MLPA Probemix P332 COL5A1 MIX-2 - 25 reactions
P332-050R	SALSA MLPA Probemix P332 COL5A1 MIX-2 - 50 reactions
P332-100R	SALSA MLPA Probemix P332 COL5A1 MIX-2 - 100 reactions
P333-025R	SALSA MLPA Probemix P333 EP300 - 25 reactions
P333-050R	SALSA MLPA Probemix P333 EP300 - 50 reactions
P333-100R	SALSA MLPA Probemix P333 EP300 - 100 reactions
P334-025R	SALSA MLPA Probemix P334 Gonadal - 25 reactions
P334-050R	SALSA MLPA Probemix P334 Gonadal - 50 reactions
P334-100R	SALSA MLPA Probemix P334 Gonadal - 100 reactions
P335-025R	SALSA MLPA Probemix P335 ALL-IKZF1 - 25 reactions
P335-050R	SALSA MLPA Probemix P335 ALL-IKZF1 - 50 reactions
P335-100R	SALSA MLPA Probemix P335 ALL-IKZF1 - 100 reactions
P336-025R	SALSA MLPA Probemix P336 UBE3A - 25 reactions
P336-050R	SALSA MLPA Probemix P336 UBE3A - 50 reactions
P336-100R	SALSA MLPA Probemix P336 UBE3A - 100 reactions
P337-025R	SALSA MLPA Probemix P337 TSC2 Confirmation - 25 reactions
P337-050R	SALSA MLPA Probemix P337 TSC2 Confirmation - 50 reactions
P337-100R	SALSA MLPA Probemix P337 TSC2 Confirmation - 100 reactions
P338-025R	SALSA MLPA Probemix P338 GBA - 25 reactions
P338-050R	SALSA MLPA Probemix P338 GBA - 50 reactions
P338-100R	SALSA MLPA Probemix P338 GBA - 100 reactions
P339-025R	SALSA MLPA Probemix P339 SHANK3 - 25 reactions
P339-050R	SALSA MLPA Probemix P339 SHANK3 - 50 reactions
P339-100R	SALSA MLPA Probemix P339 SHANK3 - 100 reactions
P340-025R	SALSA MLPA Probemix P340 EHMT1 - 25 reactions
P340-050R	SALSA MLPA Probemix P340 EHMT1 - 50 reactions
P340-100R	SALSA MLPA Probemix P340 EHMT1 - 100 reactions
P341-025R	SALSA MLPA Probemix P341 PKHD1 mix 1 - 25 reactions

P341-050R	SALSA MLPA Probemix P341 PKHD1 mix 1 - 50 reactions
P341-100R	SALSA MLPA Probemix P341 PKHD1 mix 1 - 100 reactions
P342-025R	SALSA MLPA Probemix P342 PKHD1 mix 2 - 25 reactions
P342-050R	SALSA MLPA Probemix P342 PKHD1 mix 2 - 50 reactions
P342-100R	SALSA MLPA Probemix P342 PKHD1 mix 2 - 100 reactions
P343-025R	SALSA MLPA Probemix P343 Autism-1 - 25 reactions
P343-050R	SALSA MLPA Probemix P343 Autism-1 - 50 reactions
P343-100R	SALSA MLPA Probemix P343 Autism-1 - 100 reactions
P347-025R	SALSA MLPA Probemix P347 Hemochromatosis - 25 reactions
P347-050R	SALSA MLPA Probemix P347 Hemochromatosis - 50 reactions
P347-100R	SALSA MLPA Probemix P347 Hemochromatosis - 100 reactions
P348-025R	SALSA MLPA Probemix P348 ATP1A2-CACNA1A-PRRT2 - 25 reactions
P348-050R	SALSA MLPA Probemix P348 ATP1A2-CACNA1A-PRRT2 - 50 reactions
P348-100R	SALSA MLPA Probemix P348 ATP1A2-CACNA1A-PRRT2 - 100 reactions
P350-025R	SALSA MLPA Probemix P350 CLCN1-KCNJ2 - 25 reactions
P350-050R	SALSA MLPA Probemix P350 CLCN1-KCNJ2 - 50 reactions
P350-100R	SALSA MLPA Probemix P350 CLCN1-KCNJ2 - 100 reactions
P351-025R	SALSA MLPA Probemix P351 PKD1 - 25 reactions
P351-050R	SALSA MLPA Probemix P351 PKD1 - 50 reactions
P351-100R	SALSA MLPA Probemix P351 PKD1 - 100 reactions
P352-025R	SALSA MLPA Probemix P352 PKD1-PKD2 - 25 reactions
P352-050R	SALSA MLPA Probemix P352 PKD1-PKD2 - 50 reactions
P352-100R	SALSA MLPA Probemix P352 PKD1-PKD2 - 100 reactions
P353-025R	SALSA MLPA Probemix P353 CMT4 - 25 reactions
P353-050R	SALSA MLPA Probemix P353 CMT4 - 50 reactions
P353-100R	SALSA MLPA Probemix P353 CMT4 - 100 reactions
P355-025R	SALSA MLPA Probemix P355 Microcephaly - 25 reactions
P355-050R	SALSA MLPA Probemix P355 Microcephaly - 50 reactions
P355-100R	SALSA MLPA Probemix P355 Microcephaly - 100 reactions
P357-025R	SALSA MLPA Probemix P357 MODY Mix 2 - 25 reactions
P357-050R	SALSA MLPA Probemix P357 MODY Mix 2 - 50 reactions
P357-100R	SALSA MLPA Probemix P357 MODY Mix 2 - 100 reactions
P359-025R	SALSA MLPA Probemix P359 PLOD1 - 25 reactions
P359-050R	SALSA MLPA Probemix P359 PLOD1 - 50 reactions
P359-100R	SALSA MLPA Probemix P359 PLOD1 - 100 reactions
P360-025R	SALSA MLPA Probemix P360 Y-Chromosome - 25 reactions
P360-050R	SALSA MLPA Probemix P360 Y-Chromosome - 50 reactions
P360-100R	SALSA MLPA Probemix P360 Y-Chromosome - 100 reactions
P361-025R	SALSA MLPA Probemix P361 USH2A mix 1 - 25 reactions
P361-050R	SALSA MLPA Probemix P361 USH2A mix 1 - 50 reactions
P362-025R	SALSA MLPA Probemix P362 USH2A mix 2 - 25 reactions
P362-050R	SALSA MLPA Probemix P362 USH2A mix 2 - 50 reactions
P366-025R	SALSA MLPA Probemix P366 CHM-RP2-RPGR - 25 reactions
P366-050R	SALSA MLPA Probemix P366 CHM-RP2-RPGR - 50 reactions
P366-100R	SALSA MLPA Probemix P366 CHM-RP2-RPGR - 100 reactions
P367-025R	SALSA MLPA Probemix P367 BEST1-PRPH2 - 25 reactions
P367-050R	SALSA MLPA Probemix P367 BEST1-PRPH2 - 50 reactions
P367-100R	SALSA MLPA Probemix P367 BEST1-PRPH2 - 100 reactions
P368-025R	SALSA MLPA Probemix P368 DCLRE1C - 25 reactions
P368-050R	SALSA MLPA Probemix P368 DCLRE1C - 50 reactions
P368-100R	SALSA MLPA Probemix P368 DCLRE1C - 100 reactions
P369-025R	SALSA MLPA Probemix P369 Smith-Magenis - 50 reactions
P369-050R	SALSA MLPA Probemix P369 Smith-Magenis - 50 reactions
P369-100R	SALSA MLPA Probemix P369 Smith-Magenis - 100 reactions

P370-025R	SALSA MLPA Probemix P370 BRAF-IDH1-IDH2 - 25 reactions
P370-050R	SALSA MLPA Probemix P370 BRAF-IDH1-IDH2 - 50 reactions
P370-100R	SALSA MLPA Probemix P370 BRAF-IDH1-IDH2 - 100 reactions
P377-025R	SALSA MLPA Probemix P377 Hematologic malign. - 25 reactions
P377-050R	SALSA MLPA Probemix P377 Hematologic malign. - 50 reactions
P377-100R	SALSA MLPA Probemix P377 Hematologic malign. - 100 reactions
P378-025R	SALSA MLPA Probemix P378 MUTYH - 25 reactions
P378-050R	SALSA MLPA Probemix P378 MUTYH - 50 reactions
P378-100R	SALSA MLPA Probemix P378 MUTYH - 100 reactions
P379-025R	SALSA MLPA Probemix P379 NRXN1 - 25 reactions
P379-050R	SALSA MLPA Probemix P379 NRXN1 - 50 reactions
P379-100R	SALSA MLPA Probemix P379 NRXN1 - 100 reactions
P380-025R	SALSA MLPA Probemix P380 Wilms' tumour - 25 reactions
P380-050R	SALSA MLPA Probemix P380 Wilms' tumour - 50 reactions
P380-100R	SALSA MLPA Probemix P380 Wilms' tumour - 100 reactions
P381-025R	SALSA MLPA Probemix P381 COL11A1 mix 1 - 25 reactions
P381-050R	SALSA MLPA Probemix P381 COL11A1 mix 1 - 50 reactions
P381-100R	SALSA MLPA Probemix P381 COL11A1 mix 1 - 100 reactions
P382-025R	SALSA MLPA Probemix P382 COL11A1 mix 2 - 25 reactions
P382-050R	SALSA MLPA Probemix P382 COL11A1 mix 2 - 50 reactions
P382-100R	SALSA MLPA Probemix P382 COL11A1 mix 2 - 100 reactions
P383-025R	SALSA MLPA Probemix P383 T-ALL - 25 reactions
P383-050R	SALSA MLPA Probemix P383 T-ALL - 50 reactions
P383-100R	SALSA MLPA Probemix P383 T-ALL - 100 reactions
P385-025R	SALSA MLPA Probemix P385 DOCK8 - 25 reactions
P385-050R	SALSA MLPA Probemix P385 DOCK8 - 50 reactions
P385-100R	SALSA MLPA Probemix P385 DOCK8 - 100 reactions
P386-025R	SALSA MLPA Probemix P386 DOCK8 STAT3 - 25 reactions
P386-050R	SALSA MLPA Probemix P386 DOCK8 STAT3 - 50 reactions
P386-100R	SALSA MLPA Probemix P386 DOCK8 STAT3 - 100 reactions
P387-025R	SALSA MLPA Probemix P387 NPHP1 - 25 reactions
P387-050R	SALSA MLPA Probemix P387 NPHP1 - 50 reactions
P387-100R	SALSA MLPA Probemix P387 NPHP1 - 100 reactions
P388-025R	SALSA MLPA Probemix P388 AGS - 25 reactions
P388-050R	SALSA MLPA Probemix P388 AGS - 50 reactions
P388-100R	SALSA MLPA Probemix P388 AGS - 100 reactions
P389-025R	SALSA MLPA Probemix P389 MLL2 - 25 reactions
P389-050R	SALSA MLPA Probemix P389 MLL2 - 50 reactions
P389-100R	SALSA MLPA Probemix P389 MLL2 - 100 reactions
P391-025R	SALSA MLPA Probemix P391 LAMA2 mix 1 - 25 reactions
P391-050R	SALSA MLPA Probemix P391 LAMA2 mix 1 - 50 reactions
P391-100R	SALSA MLPA Probemix P391 LAMA2 mix 1 - 100 reactions
P392-025R	SALSA MLPA Probemix P392 LAMA2 mix 2 - 25 reactions
P392-050R	SALSA MLPA Probemix P392 LAMA2 mix 2 - 50 reactions
P392-100R	SALSA MLPA Probemix P392 LAMA2 mix 2 - 100 reactions
P395-025R	SALSA MLPA Probemix P395 MEF2C-FOXP1 - 25 reactions
P395-050R	SALSA MLPA Probemix P395 MEF2C-FOXP1 - 50 reactions
P395-100R	SALSA MLPA Probemix P395 MEF2C-FOXP1 - 100 reactions
P397-025R	SALSA MLPA Probemix P397 SCN4A-CACNA1S - 25 reactions
P397-050R	SALSA MLPA Probemix P397 SCN4A-CACNA1S - 50 reactions
P397-100R	SALSA MLPA Probemix P397 SCN4A-CACNA1S - 100 reactions
P398-025R	SALSA MLPA Probemix P398 CASK - 25 reactions
P398-050R	SALSA MLPA Probemix P398 CASK - 50 reactions
P398-100R	SALSA MLPA Probemix P398 CASK - 100 reactions

P405-025R	SALSA MLPA Probemix P405 CMT1 - 25 reactions
P405-050R	SALSA MLPA Probemix P405 CMT1 - 50 reactions
P405-100R	SALSA MLPA Probemix P405 CMT1 - 100 reactions
P406-025R	SALSA MLPA Probemix P406 CMT2 - 25 reactions
P406-050R	SALSA MLPA Probemix P406 CMT2 - 50 reactions
P406-100R	SALSA MLPA Probemix P406 CMT2 - 100 reactions
P409-025R	SALSA MLPA Probemix P409 RASA1-EPHB4 - 25 reactions
P409-050R	SALSA MLPA Probemix P409 RASA1-EPHB4 - 50 reactions
P409-100R	SALSA MLPA Probemix P409 RASA1-EPHB4 - 100 reactions
P410-025R	SALSA MLPA Probemix P410 GRIN2A GRIN2B - 25 reactions
P410-050R	SALSA MLPA Probemix P410 GRIN2A GRIN2B - 50 reactions
P410-100R	SALSA MLPA Probemix P410 GRIN2A GRIN2B - 100 reactions
P411-025R	SALSA MLPA Probemix P411 Porphyria mix 1 - 25 reactions
P411-050R	SALSA MLPA Probemix P411 Porphyria mix 1 - 50 reactions
P411-100R	SALSA MLPA Probemix P411 Porphyria mix 1 - 100 reactions
P412-025R	SALSA MLPA Probemix P412 Porphyria mix 2 - 25 reactions
P412-050R	SALSA MLPA Probemix P412 Porphyria mix 2 - 50 reactions
P412-100R	SALSA MLPA Probemix P412 Porphyria mix 2 - 100 reactions
P414-025R	SALSA MLPA Probemix P414 MDS - 25 reactions
P414-050R	SALSA MLPA Probemix P414 MDS - 50 reactions
P414-100R	SALSA MLPA Probemix P414 MDS - 100 reactions
P417-025R	SALSA MLPA Probemix P417 BAP1 - 25 reactions
P417-050R	SALSA MLPA Probemix P417 BAP1 - 50 reactions
P417-100R	SALSA MLPA Probemix P417 BAP1 - 100 reactions
P418-025R	SALSA MLPA Probemix P418 MYH7 - 25 reactions
P418-050R	SALSA MLPA Probemix P418 MYH7 - 50 reactions
P418-100R	SALSA MLPA Probemix P418 MYH7 - 100 reactions
P419-025R	SALSA MLPA Probemix P419 CDKN2A/2B-CDK4 - 25 reactions
P419-050R	SALSA MLPA Probemix P419 CDKN2A/2B-CDK4 - 50 reactions
P419-100R	SALSA MLPA Probemix P419 CDKN2A/2B-CDK4 - 100 reactions
P420-025R	SALSA MLPA Probemix P420 MPN mix 1 - 25 reactions
P420-050R	SALSA MLPA Probemix P420 MPN mix 1 - 50 reactions
P420-100R	SALSA MLPA Probemix P420 MPN mix 1 - 100 reactions
P425-025R	SALSA MLPA Probemix P425 Multiple Myeloma - 25 reactions
P425-050R	SALSA MLPA Probemix P425 Multiple Myeloma - 50 reactions
P425-100R	SALSA MLPA Probemix P425 Multiple Myeloma - 100 reactions
P426-025R	SALSA MLPA Probemix P426 Cystinuria - 25 reactions
P426-050R	SALSA MLPA Probemix P426 Cystinuria - 50 reactions
P426-100R	SALSA MLPA Probemix P426 Cystinuria - 100 reactions
P429-025R	SALSA MLPA Probemix P429 SDHA-MAX-TMEM127 - 25 reactions
P429-050R	SALSA MLPA Probemix P429 SDHA-MAX-TMEM127 - 50 reactions
P429-100R	SALSA MLPA Probemix P429 SDHA-MAX-TMEM127 - 100 reactions
P431-025R	SALSA MLPA Probemix P431 FOXF1 - 25 reactions
P431-050R	SALSA MLPA Probemix P431 FOXF1 - 50 reactions
P431-100R	SALSA MLPA Probemix P431 FOXF1 - 100 reactions
P432-025R	SALSA MLPA Probemix P432 MYH9 - 25 reactions
P432-050R	SALSA MLPA Probemix P432 MYH9 - 50 reactions
P432-100R	SALSA MLPA Probemix P432 MYH9 - 100 reactions
P433-025R	SALSA MLPA Probemix P433 ARID1A-ARID1B - 25 reactions
P433-050R	SALSA MLPA Probemix P433 ARID1A-ARID1B - 50 reactions
P433-100R	SALSA MLPA Probemix P433 ARID1A-ARID1B - 100 reactions
P436-025R	SALSA MLPA Probemix P436 ANO5 - 25 reactions
P436-050R	SALSA MLPA Probemix P436 ANO5 - 50 reactions
P436-100R	SALSA MLPA Probemix P436 ANO5 - 100 reactions

P437-025R	SALSA MLPA Probemix P437 Familiar MDS-AML - 25 reactions
P437-050R	SALSA MLPA Probemix P437 Familiar MDS-AML - 50 reactions
P437-100R	SALSA MLPA Probemix P437 Familiar MDS-AML - 100 reactions
P438-025R	SALSA MLPA Probemix P438 Celiac Disease - 25 reactions
P438-050R	SALSA MLPA Probemix P438 Celiac Disease - 50 reactions
P438-100R	SALSA MLPA Probemix P438 Celiac Disease - 100 reactions
P439-025R	SALSA MLPA Probemix P439 COL4A3 - 25 reactions
P439-050R	SALSA MLPA Probemix P439 COL4A3 - 50 reactions
P439-100R	SALSA MLPA Probemix P439 COL4A3 - 100 reactions
P440-025R	SALSA MLPA Probemix P440 F10 + F11 - 25 reactions
P440-050R	SALSA MLPA Probemix P440 F10 + F11 - 50 reactions
P440-100R	SALSA MLPA Probemix P440 F10 + F11 - 100 reactions
P441-025R	SALSA MLPA Probemix P441 SACS - 25 reactions
P441-050R	SALSA MLPA Probemix P441 SACS - 50 reactions
P441-100R	SALSA MLPA Probemix P441 SACS - 100 reactions
P443-025R	SALSA MLPA Probemix P443 KANSL1 - 25 reactions
P443-050R	SALSA MLPA Probemix P443 KANSL1 - 50 reactions
P443-100R	SALSA MLPA Probemix P443 KANSL1 - 100 reactions
P444-025R	SALSA MLPA Probemix P444 COL4A4 - 25 reactions
P444-050R	SALSA MLPA Probemix P444 COL4A4 - 50 reactions
P444-100R	SALSA MLPA Probemix P444 COL4A4 - 100 reactions
P445-025R	SALSA MLPA Probemix P445 KDM6A - 25 reactions
P445-050R	SALSA MLPA Probemix P445 KDM6A - 50 reactions
P445-100R	SALSA MLPA Probemix P445 KDM6A - 100 reactions
P446-025R	SALSA MLPA Probemix P446 GALT - 25 reactions
P446-050R	SALSA MLPA Probemix P446 GALT - 50 reactions
P446-100R	SALSA MLPA Probemix P446 GALT - 100 reactions
P451-025R	SALSA MLPA Probemix P451 Chromosome 16 - 25 reactions
P451-050R	SALSA MLPA Probemix P451 Chromosome 16 - 50 reactions
P451-100R	SALSA MLPA Probemix P451 Chromosome 16 - 100 reactions
P453-025R	SALSA MLPA Probemix P453 GAA - 25 reactions
P453-050R	SALSA MLPA Probemix P453 GAA - 50 reactions
P453-100R	SALSA MLPA Probemix P453 GAA - 100 reactions
P454-025R	SALSA MLPA Probemix P454 CGD - 25 reactions
P454-050R	SALSA MLPA Probemix P454 CGD - 50 reactions
P454-100R	SALSA MLPA Probemix P454 CGD - 100 reactions
P455-025R	SALSA MLPA Probemix P455 LZTR1 - 25 reactions
P455-050R	SALSA MLPA Probemix P455 LZTR1 - 50 reactions
P455-100R	SALSA MLPA Probemix P455 LZTR1 - 100 reactions
P456-025R	SALSA MLPA Probemix P456 EVC EVC2 - 25 reactions
P456-050R	SALSA MLPA Probemix P456 EVC EVC2 - 50 reactions
P456-100R	SALSA MLPA Probemix P456 EVC EVC2 - 100 reactions
P457-025R	SALSA MLPA Probemix P457 DHCR7 - 25 reactions
P457-050R	SALSA MLPA Probemix P457 DHCR7 - 50 reactions
P457-100R	SALSA MLPA Probemix P457 DHCR7 - 100 reactions
P459-025R	SALSA MLPA Probemix P459 SERPINA1 - 25 reactions
P459-050R	SALSA MLPA Probemix P459 SERPINA1 - 50 reactions
P459-100R	SALSA MLPA Probemix P459 SERPINA1 - 100 reactions
P460-025R	SALSA MLPA Probemix P460 SMA (Silent) Carrier - 25 reactions
P460-050R	SALSA MLPA Probemix P460 SMA (Silent) Carrier - 50 reactions
P460-100R	SALSA MLPA Probemix P460 SMA (Silent) Carrier - 100 reactions
P461-025R	SALSA MLPA Probemix P461 STRC-CATSPER2-OTOA - 25 reactions
P461-050R	SALSA MLPA Probemix P461 STRC-CATSPER2-OTOA - 50 reactions
P461-100R	SALSA MLPA Probemix P461 STRC-CATSPER2-OTOA - 100 reactions

P463-025R	SALSA MLPA Probemix P463 MRKH - 25 reactions
P463-050R	SALSA MLPA Probemix P463 MRKH - 50 reactions
P463-100R	SALSA MLPA Probemix P463 MRKH - 100 reactions
P465-025R	SALSA MLPA Probemix P465 ACADM - 25 reactions
P465-050R	SALSA MLPA Probemix P465 ACADM - 50 reactions
P465-100R	SALSA MLPA Probemix P465 ACADM - 100 reactions
P466-025R	SALSA MLPA Probemix P466 CDC73 - 25 reactions
P466-050R	SALSA MLPA Probemix P466 CDC73 - 50 reactions
P466-100R	SALSA MLPA Probemix P466 CDC73 - 100 reactions
P469-025R	SALSA MLPA Probemix P469 F5 - 25 reactions
P469-050R	SALSA MLPA Probemix P469 F5 - 50 reactions
P469-100R	SALSA MLPA Probemix P469 F5 - 100 reactions
P470-025R	SALSA MLPA Probemix P470 NCL - 25 reactions
P470-050R	SALSA MLPA Probemix P470 NCL - 50 reactions
P470-100R	SALSA MLPA Probemix P470 NCL - 100 reactions
P471-025R	SALSA MLPA Probemix P471 EOFAD - 25 reactions
P471-050R	SALSA MLPA Probemix P471 EOFAD - 50 reactions
P471-100R	SALSA MLPA Probemix P471 EOFAD - 100 reactions
P472-025R	SALSA MLPA Probemix P472 SUFU - 25 reactions
P472-050R	SALSA MLPA Probemix P472 SUFU - 50 reactions
P472-100R	SALSA MLPA Probemix P472 SUFU - 100 reactions
P473-025R	SALSA MLPA Probemix P473 CTNS - 25 reactions
P473-050R	SALSA MLPA Probemix P473 CTNS - 50 reactions
P473-100R	SALSA MLPA Probemix P473 CTNS - 100 reactions
P474-025R	SALSA MLPA Probemix P474 CD274-PDCD1LG2-JAK2 - 25 reactions
P474-050R	SALSA MLPA Probemix P474 CD274-PDCD1LG2-JAK2 - 50 reactions
P474-100R	SALSA MLPA Probemix P474 CD274-PDCD1LG2-JAK2 - 100 reactions
P475-025R	SALSA MLPA Probemix P475 FOXP1 - FOXP2 - 25 reactions
P475-050R	SALSA MLPA Probemix P475 FOXP1 - FOXP2 - 50 reactions
P475-100R	SALSA MLPA Probemix P475 FOXP1 - FOXP2 - 100 reactions
P476-025R	SALSA MLPA Probemix P476 ZNRF3 - 25 reactions
P476-050R	SALSA MLPA Probemix P476 ZNRF3 - 50 reactions
P476-100R	SALSA MLPA Probemix P476 ZNRF3 - 100 reactions
P478-025R	SALSA MLPA Probemix P478 SMARCE1 - 25 reactions
P478-050R	SALSA MLPA Probemix P478 SMARCE1 - 50 reactions
P478-100R	SALSA MLPA Probemix P478 SMARCE1 - 100 reactions
P479-025R	SALSA MLPA Probemix P479 TCF12 - ERF - 25 reactions
P479-050R	SALSA MLPA Probemix P479 TCF12 - ERF - 50 reactions
P479-100R	SALSA MLPA Probemix P479 TCF12 - ERF - 100 reactions
P480-025R	SALSA MLPA Probemix P480 WHS & Achondroplasia - 25 reactions
P480-050R	SALSA MLPA Probemix P480 WHS & Achondroplasia - 50 reactions
P480-100R	SALSA MLPA Probemix P480 WHS & Achondroplasia - 100 reactions
P481-025R	SALSA MLPA Probemix P481 PRKAR1A-ARMC5 - 25 reactions
P481-050R	SALSA MLPA Probemix P481 PRKAR1A-ARMC5 - 50 reactions
P481-100R	SALSA MLPA Probemix P481 PRKAR1A-ARMC5 - 100 reactions
P482-025R	SALSA MLPA Probemix P482 DICER1 - 25 reactions
P482-050R	SALSA MLPA Probemix P482 DICER1 - 50 reactions
P482-100R	SALSA MLPA Probemix P482 DICER1 - 100 reactions
P483-025R	SALSA MLPA Probemix P483 HER gene family - 25 reactions
P483-050R	SALSA MLPA Probemix P483 HER gene family - 50 reactions
P483-100R	SALSA MLPA Probemix P483 HER gene family - 100 reactions
P484-025R	SALSA MLPA Probemix P484 ALPL - 25 reactions
P484-050R	SALSA MLPA Probemix P484 ALPL - 50 reactions
P484-100R	SALSA MLPA Probemix P484 ALPL - 100 reactions

P488-025R	SALSA MLPA Probemix P488 RS1 - 25 reactions
P488-050R	SALSA MLPA Probemix P488 RS1 - 50 reactions
P488-100R	SALSA MLPA Probemix P488 RS1 - 100 reactions
P489-025R	SALSA MLPA Probemix P489 BARD1 - 25 reactions
P489-050R	SALSA MLPA Probemix P489 BARD1 - 50 reactions
P489-100R	SALSA MLPA Probemix P489 BARD1 - 100 reactions
P490-025R	SALSA MLPA Probemix P490 ADA2 - 25 reactions
P490-050R	SALSA MLPA Probemix P490 ADA2 - 50 reactions
P490-100R	SALSA MLPA Probemix P490 ADA2 - 100 reactions
P492-025R	SALSA MLPA Probemix P492 POLD1 - POLE - 25 reactions
P492-050R	SALSA MLPA Probemix P492 POLD1 - POLE - 50 reactions
P492-100R	SALSA MLPA Probemix P492 POLD1 - POLE - 100 reactions
p493-025R	SALSA MLPA Probemix p493 SMCHD1 - 25 reactions
p493-050R	SALSA MLPA Probemix p493 SMCHD1 - 50 reactions
p493-100R	SALSA MLPA Probemix p493 SMCHD1 - 100 reactions
P494-025R	SALSA MLPA Probemix P494 NBN - 25 reactions
P494-050R	SALSA MLPA Probemix P494 NBN - 50 reactions
P494-100R	SALSA MLPA Probemix P494 NBN - 100 reactions
P495-025R	SALSA MLPA Probemix P495 CYP11A1-CYP11B1-CYP11B2 - 25 reactions
P495-050R	SALSA MLPA Probemix P495 CYP11A1-CYP11B1-CYP11B2 - 50 reactions
P495-100R	SALSA MLPA Probemix P495 CYP11A1-CYP11B1-CYP11B2 - 100 reactions
P520-025R	SALSA MLPA Probemix P520 MPN mix 2 - 25 reactions
P520-050R	SALSA MLPA Probemix P520 MPN mix 2 - 50 reactions
P520-100R	SALSA MLPA Probemix P520 MPN mix 2 - 100 reactions